



TopLine Sanitisers

Active protection of patients
and nursing staff






DISTRIBUTED BY

What do people all over the world expect from today's nursing profession? Certainly compliance with relevant laws and standards and a duty of care, but also proper adherence to occupational health and safety guidelines plus regular monitoring and documentation of infection control measures. The aim is clear – to prevent infectious diseases in hospitals and care homes.

MEIKO can actively assist you in this task by providing professional expertise and reliable cleaning and disinfection technology, so you get the reassurance of reliable hygiene from a single source. We can help you plan your facilities and choose the right technology and chemicals. What's more, we offer superb technical support and professional consulting services even after the sale, because we believe in forging long-term partnerships with our customers. Whatever your needs, our bespoke hygiene management solutions have got you covered. Meeting statutory and infection control requirements is a serious business – but with MEIKO at your side you can be confident of getting it right.

Here's how MEIKO TopLine helps prevent infections and protect patients and staff:

- Disinfects the entire system every time you start a wash program
- Cools and dries care utensils with HEPA-filtered air
- Completely steam-tight system. No steam is released, not even when the door is opened at the end of the wash program
- Certified to DIN EN ISO 15883 by an independent test institute which, in turn, is accredited by the central authority of the German federal Länder for health protection regarding medicinal products and medical devices
- Energy-optimised disinfection process. The machine only uses the energy required to reach the preset A_0 value
- Displays target and actual A_0 value
- Allows user to select A_0 value (A_0 60 / A_0 600 / A_0 3000)
- Helpful and informative 4 x 20-character plain text display with energy-saving mode
- Nozzle system configured to get care utensils perfectly clean, including a unique controllable telescopic rotary jet designed to thoroughly clean the inside of the utensils



Optimises hygiene and safety with a patented design which ensures that no steam escapes from the system.

Gets every last spot of your care utensils hygienically clean with a powerful telescopic rotary jet and an optimised wash system that acts on all sides.

Takes reliable disinfection to the next level by offering patented control of the freely selectable* A_0 value, carefully monitored by temperature sensors which cross-check each other's readings.

*from 60 to 3000.

Uses filtered air to dry and cool, producing hygienically clean care utensils that are cool to the touch.

The *Total Auto Disinfection (T.A.D.)* function disinfects all the water pipes in the machine after each completed wash cycle to avoid the risk of re-contamination.

Uses a unique suction lance coding system to eliminate the risk of mixing up rinse aid and detergent and to ensure a validated process.

Guarantees high operating reliability by using sensors to monitor dosing of the detergent and of the softener/rinse aid mix.

Many of the dirty utility rooms of hospitals and care homes no longer comply with current regulations and hygiene standards as well as the demands of nursing staff. TopLine single appliances were specially designed for rapid and uncomplicated installation to replace out-dated cleaning and disinfection technology.

TopLine 10 & TopLine 20 - for the clear-of-the-floor wall mounting & floor mounted appliances



TopLine 10

- Clear-of-the-floor wall mounting
- The easy way of installing a cleaning and disinfection appliance.
- Space-saving and ready for connection.
- A practical and compact work unit
- Available with automatic doors



TopLine 20

- Floor mounted appliances
- Cleaning and disinfection appliance with plinth. Just place in position and connect.
- The ideal way of utilising existing connections.
- A work unit which is versatile in its functions, easy to install and compact.
- Available with automatic doors

Code	Description	Width	Depth	Height
UPM200100	Sanitiser - Meiko - TopLine 10	500 mm	450 or 588 mm	1630 mm
UPM200110	Sanitiser - Meiko - TopLine 10 AT - with Auto Door	500 mm	450 or 588 mm	1630 mm
UPM200120	Sanitiser - Meiko - TopLine 20	500 mm	450 - 588 mm	1730 mm
UPM200130	Sanitiser - Meiko - TopLine 20 AT - with Auto Door & Foot	500 mm	450 - 588 mm	1730 mm

They can be modified to fit into the position of the old appliance whatever the make of the original - and to use the same connections.

TopLine 20XL - floor mounted appliances



TopLine 20XL

- Double work unit which is versatile in its functions, easy to install and compact.
- Floor mounted appliances
- Cleaning and disinfection appliance with plinth. Just place in position and connect.
- The ideal way of utilising existing connections.
- Available with automatic doors

Code	Description	Width	Depth	Height
UPM200140	Sanitiser - Meiko - TopLine 20 XL (side-by-side)	1000 mm	450 - 588 mm	1730 mm
UPM200150	Sanitiser - Meiko - TopLine 20 XL AT (side-by-side) - with Auto Door & Foot	1000 mm	450 - 588 mm	1730 mm

MEIKO TopLine single appliances – seamless and adaptable hygiene security for every work room

Accessories and range of fixtures



TopLine 30 & TopLine 30 with WC Outlet - for the flush wall installation into bathrooms



TopLine 30

- Single wall mounted appliance – seamless and adaptable
- Hygiene security for every work room
- Cleaning and disinfection appliance TopLine 30 for the bathroom of the patients' rooms, e.g. on the isolation ward.
- With mounting frame
- Available with automatic doors

Code	Description	Basic Frame	Visible Fronts
UPM200160	Sanitiser - Meiko - TopLine 30	650 x 450 x 2210 mm	530 x 1760 mm
UPM200170	Sanitiser - Meiko - TopLine 30 AT - with Auto Door	650 x 450 x 2210 mm	530 x 1760 mm

TopLine 30 with WC Outlet

- Cleaning and disinfection appliance TopLine 30 with WC outlet for the bathroom in the patients' room, e.g. on the isolation ward.
- With mounting frame and complete piping as well as central connections for all commercially available WC types.
- Available with automatic doors

Code	Description	Basic Frame	Visible Fronts
UPM200180	Sanitiser - Meiko - TopLine 30 WC	650 x 450 x 2210 mm	530 x 1760 mm
UPM200190	Sanitiser - Meiko - TopLine 30 WCAT - with Auto Door	650 x 450 x 2210 mm	530 x 1760 mm

In addition to the Meiko cleaning and disinfection technology, we can also supply you with a comprehensive range of accessories for safe and ergonomic working and the hygienic storage of care utensils.



Accessories

	Code	Description	Width	Depth	Height
1	UPA400100	Basket insert with hinged lid - Meiko	-	-	-
2	UPA400110	Basket with hinged lid for small items - Meiko	-	-	-
3	UPA400200	Bottle Rack - Meiko TU50 - 8 Urine Bottles (drip tray sep)	500 mm	350 mm	650 mm
4	UPA400210	Bottle Rack - Meiko TU70 - 11 Urine Bottles (drip tray sep)	700 mm	350 mm	650 mm
5	UPA400250	Pan & Bottle Rack - Meiko TUS 140 - 9 Pan & 9 Bottle (drip tray sep)	1400 mm	350 mm	650 mm
6	UPA400260	Pan & Bottle Rack - Meiko TUS 50 - 6 Pan & 6 Bottle (drip tray sep)	500 mm	350 mm	650 mm
7	UPA400300	Pan Rack - Meiko TS100 - 12 Pans (drip tray sep)	1000 mm	350 mm	650 mm
8	UPA400310	Pan Rack - Meiko TS50 - 3 Pans (drip tray sep)	500 mm	350 mm	650 mm
9	UPA400320	Pan Rack - Meiko TS70 - 6 Pans (drip tray sep)	700 mm	350 mm	650 mm
10	UPA400330	Pan Rack - Meiko TS90 - 9 Pans (drip tray sep)	900 mm	350 mm	650 mm
11	UPA400340	Pan Rack (Not Mounted) - Meiko - For vertical storage of bedpans - 4 Pan	-	-	-
12	UPA400350	Pan Rack (Not Mounted) - Meiko - For vertical storage of bedpans - 8 Pan	-	-	-
13	UPA400400	Rack (Internal) - Meiko - SUH 1 Internal rack (Non-standard in Aust)	-	-	-
14	UPA400410	Rack (Internal) - Meiko - SUH 2 Quick-change rack (Standard in Aust)	-	-	-
15	UPA400420	Rack (Internal) - Meiko - TopCut B - Quick Change Rack with Bag Empty	-	-	-
16	UPA400430	Rack (Internal) - Meiko - TopCut H - Bag Emptying Rack (internal)	-	-	-
17	UPA400500	Drip tray (For all wall-mounted racks) - Meiko - Stainless Steel (order with rack)	-	-	-



ABR95.3 190721



The clean solution

DISTRIBUTED BY



1300 133 120

www.aidacare.com.au