



# MAXI BEDSIDE COMMUNE

Experience enhanced convenience and comfort with the Aspire MAXI bedside commode. Crafted for bariatric users with a generously proportioned seat and extensive width adjustment between the arms, it provides ample space for larger users up to 310kg SWL. The soft padded seat cushion, toilet seat, backrest and armrests provide both comfort and support when sitting down, toileting and standing back up, particularly with prolonged use.



## FEATURES & BENEFITS

- **Promoting Independence:** Helps to maintain the users independence within their own home by minimizing risky night-time walks to the bathroom and the need for carer assistance, promoting independence and self-care.
- **Superior Durability:** Constructed with high-quality materials, it can withstand substantial weight, making it a reliable choice for MAXI users. The robust frame ensures stability, providing users with a sense of security and peace of mind.
- **Ease of Use:** Soft padded seat is easily lifted for toileting. Flip-up toilet seat for bowl access and cleaning.
- **Adjustable and Versatile:** Height adjustable legs and width adjustable frame means it is easily adjusted to suit a wide range of user's sizes and needs.
- **Bedroom Convenience:** Toileting and a comfortable chair are both conveniently and safely accessible within the bedroom.



CHAMPAGNE VINYL



BEDSIDE  
COMMUNE  
BEC045125

OVERALL HEIGHT mm/in	OVERALL WIDTH mm/in	OVERALL DEPTH mm/in	SEAT HEIGHT mm/in	SEAT WIDTH mm/in	SEAT DEPTH mm/in	BETWEEN ARMS mm/in	BACKREST HEIGHT mm/in	UNIT WEIGHT kg/lbs
890-1015 / 35-40	715-915 / 28-36	645 / 25.5	490-615 / 19.5-24	530 / 21	475 / 18.5	575-775 / 22.5-30.5	400 / 16	12.3 / 27

ASFY108.7180.424