



SOCIALITE WHEELCHAIR

Arrive in style

The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends. The generous 115kg safe working load caters to a broad range of users.

Padded and breathable seat and back provide optimum comfort and postural support. Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.

LESS THAN
9kg /
19.8lbs



FEATURES & BENEFITS

- Sleek European design
- Stylish contemporary colours
- One of the lightest transport chairs available
- Perfect for outings and appointments
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Contoured seating improves comfort and postural support
- Rear storage pocket for essentials
- Available in attendant (12 in) and self-propelled (22 in) options



Attendant Propelled

Self Propelled



UPHOLSTERY
Red/Black
FRAME
Red



UPHOLSTERY
Silver/Black
FRAME
Silver

1 YEAR WARRANTY
Components

5 YEAR WARRANTY
Frame



ALUMINIUM

	OVERALL WIDTH (mm/inches)	SEAT DEPTH (mm/inches)	SEAT WIDTH (mm/inches)	SEAT HEIGHT ANTERIOR (mm/inches)	SEAT HEIGHT POSTERIOR (mm/inches)	BACKREST HEIGHT (mm/inches)	ARMREST HEIGHT (mm/inches)	PRODUCT WEIGHT (kg / lbs)	MIN. CHAIR WEIGHT (kg / lbs)
Attendant Propelled MWS449800 - Silver MWS449810 - Red	675 / 26.6	405 / 16	450 / 17.7	475-500 / 18.7-19.7	450 / 17.7	410 / 16.1	200 / 7.9	9 / 19.8	8 / 17.6
Self Propelled MWS449820 - Silver MWS449830 - Red	675 / 26.6	405 / 16	450 / 17.7	475-500 / 18.7-19.7	450 / 17.7	410 / 16.1	200 / 7.9	13 / 28.7	8 / 17.6

ASFY7616 090224