




BEDROOM SAFETY & SUPPORT

The Aspire bedroom safety & support range has been designed with user support and safety in mind. Bed poles, cradles and rails ensure carers can provide safety and support for residents without compromising comfort in their bedroom environment.




Drop Down Bed Rail
BEA013700



Bed Pole With Adjustable
Crook Handle
BEA014150



Bed Pole With Adjustable
Crook Handle
BEA013100



Fixed Bed Rail
BEA013600



Clamp On Bed Pole
BEA009300 & BEA009310



Adjustable Height Bed Cradle
BEA014610



Fixed Height Bed Cradle
BEA014605

✓ DVA CONTRACTED

✓ Bedroom Accessories are registered
in the ARTG, AUST R ARTG 290769

1 YEAR WARRANTY



FREE STANDING SELF HELP POLE

BEA013100

- Free standing support handle provides assistance for the user to reposition or transfer in and out of bed
- Handle height is easily adjusted to suit the user
- Sits behind or beside the bed



BED POLE WITH ADJUSTABLE CROOK HANDLE

BEA014150

- Support pole designed to assist the user to balance while getting into and out of bed
- Adjustable crook handle can be located in multiple positions
- Base is positioned under the mattress



CLAMP ON BED POLE

BEA009300 (RIGHT) & BEA009310 (LEFT)

- Support pole clamps on to most bed frames to provide assistance with repositioning or transferring in and out of bed
- Folding crook handle allows the pole to swing down and out of the way as required
- Clamps onto the side of most bed frames



FIXED BED RAIL

BEA013600

- Support and safety rail provides repositioning and transfer support and assists to prevent falls from bed
- Base is positioned under the mattress



DROP DOWN BED RAIL

BEA013700

- Support and safety rail when raised is effective at preventing falls from bed
- May be lowered down and out of the way for patient transfers or repositioning
- Base is positioned under the mattress



FIXED HEIGHT BED CRADLE

BEA014605

- Keeps bed linen off sensitive areas of the patient by creating clearance between linen and mattress
- Base is positioned under the mattress



ADJUSTABLE HEIGHT BED CRADLE

BEA014610

- Keeps bed linen off sensitive areas of the patient by creating clearance between linen and mattress
- Adjustable height suitable for a range of mattress heights
- Base is positioned under the mattress