

PRODUCT UPDATE



Aspire Altitude Chair Range – Power Cord Update

We have updated the length of the power cord on the Altitude Chair range from 1.5 to 2.15 metres.

This change ensures customers can have their Altitude Chair positioned sufficiently away from the wall whilst still plugged into mains power and charging.

Additionally, we have added an upholstery label to the rear of the chair, which serves as a reminder to customers to maintain a minimum distance of 550mm between the chair and wall.



CHP227700
Altitude Vertical Lift Chair



CHP227705
Altitude Recliner Chair