

PRIDE MOBILITY SCOOTERS - CHOICE ASSISTANCE GUIDE

SCOOTER RANGE OVERVIEW



Model	Boot Scooter	Boot Scooter Supalite	Go Go LX + CTS	Sonic Scooter	Quest Folding Compact	Gopher DX Sport	GT Scooter	Pathrider 10	Pathrider 10 DX
Size	*SMALL*	*SMALL*	*SMALL*	*SMALL*	*SMALL*	**MEDIUM**	**MEDIUM**	**MEDIUM**	**MEDIUM**
Code	SCT670613	SCT670616	SCT670612	SCT670665	SCT670630	SCT670638	SCT670670	SCT670615	SCT670620
Range (Max.)	11.5 km	14.4 km	17 km	20 km	20 km	33 km	30 km	40 km	33 km
Battery	2x 12V 12AH	1x 24V 12AH	2x 12V 18AH	2x 12V 20AH	1x 24V 11.5AH Li-Polymer	2x 12V 50AH	2 x 12V 35AH	2x 12V 31AH (40AH optional)	2x 12V 31AH (40AH optional)
Seat Type	Low back	Low back	Low back	Low Back	Low back	High back	Low Back	Low back	High back
Max Weight	118 kg	120 kg	136 kg	136 kg	115 kg	160 kg	159 kg	159 kg	159 kg
Std Body Colours	Red	White	Red	Red	White	Blue	Red	Viper Blue	Viper Blue
Features	Lightweight, disassembles , dual-voltage charger, adjustable tiller, flat free tyres, front basket	Lightweight Folding design, disassembles into 6 lightweight pieces	Lightweight, disassembles , comfort track suspension (CTS), adjustable tiller, front basket	Lightweight, disassembles , LED light, adjustable tiller and delta handle	Lightweight Folding design, folds in three easy steps, adjustable tiller	High performance, sporty, flat free tyres, head & armrests, front & rear suspension	Advanced suspension, LED light, deluxe tiller, directional signals	High performance, sporty, flat free tyres, back lit battery gauge, armrests, front basket	High performance, sporty, flat free tyres, head & armrests, front & rear suspension

SCOOTER RANGE OVERVIEW



Model	Victory 10 LX Scooter	Pursuit 130 XL	Pathrider 140 XL	Celebrity Deluxe XL	Sportrider Scooter	Rangerider Scooter
Size	**MEDIUM**	***LARGE***	***LARGE***	***LARGE***	***LARGE***	***LARGE***
Code	SCT670640	SCT670571	SCT670576	SCT670540	SCT670585	SCT670581
Range (Max.)	26 km	25 km	35 km	40 km	50 km	36 km
Battery	2x 12V 40AH	2x 12V 55AH	2x 12V 70AH	2x 12V 40AH	2x 12V 70AH	2x 12V 100AH
Seat Type	High back	CRS Deluxe	Deluxe Captains Seat	Deluxe reclining high back	High back, adjustable	Pillow-top, high back
Max Weight	181 kg	159 kg	181 kg	159 kg	181 kg	159 kg
Std Body Colours	Red	Candy Apple Red	Candy Apple Red, Silver & Yellow	Candy Apple Red & Silver	Silver	Stone Grey
Features	Sporty, outdoor use, digital display speedometer/odometer. Full suspension, front/rear lighting kit	High performance, sporty, large pneumatic tyres, head & armrests, full lighting package	Top of range sporty, large pneumatic tyres for outdoor, head & armrests, front & rear suspension	High performance, large tyres for outdoor, head & armrests, lumbar support, rear suspension	Sleek outdoor scooter, digital speedometer display, full suspension system, front/rear lighting	Adjustable suspension, dual hydraulic brakes, built in USB charger, full lighting package

SCOOTER ACCESSORIES

- A. Rear Scooter Bag
- B. Scooter Flag
- C. Sun Canopy
- D. Walking Frame Holder



All information within this document is correct at the time of publication. This document should be used as a guide only. For further enquiries regarding product specifications, please refer to product brochures.
 *Range Per Charge varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition & tyre condition. Battery upgrades (where possible) improve Range Per Charge.
 **All prices are recommended retail. This varies with promotions. Contact Aidacare for the most up to date pricing.