


# Premium Pads & Mats



## FEATURES

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura<sup>1</sup> Bed/Chair Alert monitors
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress
- Floors mats made from heavy duty non-slip material

## TECHNICAL SPECIFICATIONS

- Premium Bed / Chair Pads**
- 2626 Premium Bed Pad (760x510) 1 Year
  - 2627 Premium Bed Pad (760x255) 1 Year
  - 2634 Premium Bed Pad (760x255) 30 Day
  - 2628 Premium Chair Pad (380x255) 1 Year
  - 2635 Premium Chair Pad (380x255) 30 Day
- Premium Toilet Pad**
- 2631 Premium Toilet Pressure Pad (140x40) 30 Day
- Premium Floor Mats**
- 2633 Premium Floor Mat (1220x610) 1 Year
  - 2637 Premium Floor Mat (1625X610) 1 Year
  - Weight range for reliable operation: 30kg - 200kg
  - Bottom material PVC without Phthalate.

## ACCESSORIES

**2611** Standard Monitor, **2607** Safety Reset Monitor, **2608** Wireless Monitor, **2613** Voice Monitor, **2614** SafeTime Monitor\*

\*WHEN USED WITH PREMIUM MATS AND PADS THE SAFE TIME ALERT MONITOR WILL NOT HAVE SAFE TIME ALERT FEATURES