

# Q6 EDGE 3 *Stretto*<sup>™</sup>



Q6 EDGE 3 STRETTO<sup>™</sup>  
shown with iLevel<sup>®</sup> power seat elevation,  
Q-Logic 3 environmental controls, and  
Stealth<sup>®</sup> TRU-Comfort 2 cushion.



The Q6 Edge 3 Stretto<sup>™</sup> features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 31.8cm (12.5") or 35.6cm (14") drive wheels and optional 30.5cm (12") of iLevel<sup>®</sup> power adjustable seat height.



Narrowest  
complex  
rehab base  
on the market

  
**QUANTUM<sup>®</sup>**  
#1 FOR REHAB POWER



[quantumrehab.com](http://quantumrehab.com)

## Specifications

<b>DRIVE WHEELS</b>	35.6cm (14") or 31.8cm (12.5") pneumatic or flat free tires
<b>CASTER WHEELS</b>	<b>Front:</b> 12.7cm (5") flat free <b>Rear:</b> 15.2cm (6") flat free
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
<b>GROUND CLEARANCE</b>	7.4cm (2.9")
<b>TURNING RADIUS</b> <sup>2</sup>	53.3cm (21")
<b>OVERALL LENGTH</b> <sup>2</sup>	82.5cm (32.5") without front riggings
<b>BASE WIDTH</b>	55.25cm (21.8") with 35.6cm (14") drive wheels 52.71cm (20.8") with 31.8cm (12.5") drive wheels

### SEATING SIZES

<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12"W-16"W <b>Depth:</b> 12"W-16"W
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12"W-16"W <b>Depth:</b> 12"W-16"W
<b>TRU-Balance® 3 Tilt and iLevel®</b>	<b>Width:</b> 12"W-15"W OR 15"W-16"W <b>Depth:</b> 12"D-16"D OR 16"D-20"D

### SEAT-TO-FLOOR HEIGHTS

<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	40.6cm (16")-45.7cm (18")
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	44.5cm (17.5")-49.5cm (19.5")
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	44.5cm (17.5")-49.5cm (19.5")

<b>MANUAL RECLINE</b> <sup>6</sup>	Yes (TRU-Balance® Flex only)
<b>BATTERY SIZE</b> <sup>3</sup>	40 AH (2) or 55 AH (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	40 AH: 14 kg each 55 AH: 18.1 kg each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 EX
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole, Mid-Wheel 6®
<b>WEIGHT CAPACITY</b> <sup>4</sup>	136 kg
<b>BASE WEIGHT WITHOUT SEATING SYSTEM</b>	57.3 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 0.3cm (0.1") with 45.7cm (18") x 45.7cm (18") seat, fixed towers, 136 kg weight capacity and solid tires (add 0.95cm (0.4") with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 1.4 kg based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.

## Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5cm (12") of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included
- Assembled in Australia of foreign and domestic parts
- WC 19 Occupied transit approved

## Options

- iLevel® option allows 30.5cm (12") of seat elevation while at walking speed (5.6 km/h)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- 30.5cm (12") x 30.5cm (12") seat sizes with iLevel (SOF only)

## Colours\*

