

## **MOTOMed Muvi**





The unique MOTOmed muvi feature enables leg and arm training at the same time. For the first time it is possible to train both areas simultaneously. The passive or active training, movement speed and the resistance level for legs and arms can be adjusted separately for simultaneous movement. Consequently, the legs can be engaged in passive training while the user moves actively with the arms at the same time. Depending on the indication, simultaneous movement can improve the therapy effect and reduces the therapy time.

- Passive and active training modes
- Allows for simultaneous training of legs and arms in both passive or active modes
- Serial port activated: trainer can be upgraded to an FES Cycle with the HASOMED RehaMove2 'Functional Electrical Stimulator'
- Electronic foot insertion aid
- Speed from 1 to 60 rpm (passive)
- Resistance from 0 to 20 (active)
- Motor power from 1 to 10 (passive)
- Presetting of therapy time from 0 to 120 minutes
- Movement Protection & Spasm Control
- Safety stop
- SymmetryTraining of leg and arm/upper body trainer
- Quiet, gentle and harmonious cycle movement (SmoothDriveSystem)
- Gentle starting and stopping of the movement
- Therapy and motivation programs
- Detailed biofeedback during and after training
- Export of training data and import of own training programs
- Games & videos
- Slideshow via USB (delivered without USB stick)





Language selection

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101100 N	AOTOMed Muvi - Leg guides. Foot & Hand Fix	95 - 108 cm	60 cm	124 - 135 cm	135 kg	51 kg