


Invacare Pivot Sling



	Shoulder Width (mm)	COCCYX TO TOP OF HEAD (mm)
S	600 - 700	640 - 780
M	700 - 810	750 - 900
L	760 - 870	830 - 1000



- The sling allows a patient to be lifted in either a sitting or lying position
- The ability to position a patient in the sitting posture greatly assists correct seating
- A more open lift. Recommended for general lifting including floor lifting

Warranty Information

- 1 Year



ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS374600	Invacare Pivot Mesh - Small	205 kg	Mesh
LSS374700	Invacare Pivot Mesh - Medium	205 kg	Mesh