

## Quickie M6 Bariatric Wheelchair



- Modular, High-Strength Frame
- Locking Crossbrace and Stabilizer Bar
- Heavy duty components accommodate changing conditions with simple parts replacements
- Adjustable axle plate allows for better wheel access, while low seat heights increase access to the floor for foot propulsion



Code	Display Name	Total Width	Safe Working Load	Unit Weight	Seat Depth	Seat Width
MWC427325	Quickie M6 Manual Wheelchair (As per script)	690 - 980 mm	295 kg	24 kg	450 -550 mm	560 - 750 mm