

## Easy Slide Donning Aid

- Application aids for both open-toe and closed-toe stockings
- Also suitable for hospital stockings with inspection hole
- Slide stocking over, and then remove through open toe
- Suggest to use with rubber gloves



### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### ndis

#### Item Code:

05\_090900046\_0103\_1\_2

Click [here](#) for more information



### Compliance

No applicable compliance standard

Code	Display Name
GMS340952	EasySlide for O/Toe Stockings - Med
GMS340953	EasySlide for O/Toe Stockings - Lge
GMS340954	EasySlide for O/Toe Stockings - XL