

## Batteries & Chargers



- A range of durable replacement batteries for Mobility Scooters and Wheelchairs
- Please check your manufacturer's documentation for the correct battery

- ✔ Contracted DVA product
- ✔ Contract approved

### Warranty Information

- 6 Months

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height
MBT404100	Lithium Ion Battery for HS268 - 11.6 Amp/Hour	999	1	-	-	-
MBT410000	Sealed Gel Battery - 70/73 Amp/Hour	729	1	275 mm	170 mm	235 mm
MBT410005	Sealed Gel Battery - 50 Amp/Hour - for Power Wheelchairs only	509	1	245 mm	138 mn	235 mm
MBT410007	Sealed Lead Acid Battery - 50 Amp/Hour - for Mobility Scooters only	489	1	199 mm	166 mm	171 mm
MBT410010	Sealed Gel Battery - 60 Amp/Hour	659	1	260 mm	170 mm	180 mm
MBT410015	Sealed Gel Battery - 40 Amp/Hour	463	1	197 mm	167 mm	175 mm
MBT410025	Sealed Gel Battery - 12 Amp/Hour	133	1	151 mm	97 mm	100 mm

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height
MBT410030	Sealed Gel Battery - 33 Amp/Hour	371	1	197 mm	167 mm	175 mm
MBT410035	Sealed Gel Battery - 21 Amp/Hour	212	1	180 mm	75 mm	167 mm
MBT410040	Sealed Gel Battery - 80 Amp/Hour	609	1	320 mm	170 mm	230 mm
MBT410045	Battery - 100 Amp/Hour - MK Sealed Lead Acid – for Mobility Scooters	919	1	320 mm	170 mm	235 mm
MBT403100	Lithium-Ion Battery for Pride Quest Scooter - 24V	-	-	-	-	-