

QLD HEALTH EQUIPMENT CATALOGUE

JULY 2019

VERSION: 1 UPDATED: JULY 2019



Quality Equipment Enterprise Solutions Home & Community





Aidacare & QLD Health

Aidacare are contracted to QLD Health for the supply of a range of beds and mattresses suitable for use in hospitals and health services within QLD. Aidacare have been a reliable and trusted partner for similar contracts for more than 20 years, providing National and state-wide coverage for the supply of equipment and services.

Our electronic catalogue is available for download at

www.aidacare.com.au/globalassets/resources

For information on any of the extensive range of products available from Aidacare, please browse the complete online Aidacare catalogue at www.aidacare.com.au

For additional information, please contact your local Aidacare office directly on 1300 133 120 or email salesbris@aidacare.com.au

Major Cities

State Wide Distribution Network

Aidacare



BRISBANE SOUTH

30 Dulacca Street Acacia Ridge QLD, 4110 Tel: 07 3086 2900 Fax: 07 3345 7961



BRISBANE NORTH

354 Bilsen Road Geebung QLD, 4034 Tel: 07 3086 2950 Fax: 07 3865 7896



GOLD COAST

Units 1 & 2, 48 Olympic Circuit, Southport QLD, 4215 Tel: 07 5532 6070 Fax: 07 5503 1803



TOOWOOMBA

3/342 Taylor Street Toowoomba QLD, 4350 Tel: 07 3086 2990 Fax: 07 4632 6473



SUNSHINE COAST

Unit 2, 35 Premier Circuit Warana, QLD, 4575 Tel: 07 3086 2980 Fax: 07 3345 7961



HERVEY BAY

Unit 3, 19 Islander Road, Hervey Bay, QLD, 4655 Tel: 07 4194 5019



CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).

Table of Contents



Product Index

| Home Care & Hospital Beds | Pg 04 | (1) |
|----------------------------------|-------|-----|
| Aidacare FL250 Floorline Bed | | |
| Mattress Overlays & Replacements | Pg 06 | (1) |
| Malvestio | | |







Bed - Aidacare FL250 Floorline Bed - Single - 4 Function - 250kg



- True floor bed with a low height of 110mm
- Brake First Automatic castor lock at the bed's lowest height
- Exit Safe One touch button activates CPR & safe handling bed light
- Backlit handset with seating and CPR function
- Wall buffer/protection bar
- Attendant control panel
- Bed extension for taller residents
- Mattress retainers
- Auto-regression backrest function
- Underbed lighting
- 250kg Safe Working Load

| Code |
|---------------------------------|
| Overall Length |
| Overall Width |
| Height Range |
| Backrest Incline / Decline |
| Knee Bend Incline / Decline |
| Underbed Hoist Clearance |
| SWL |
| |

| ranty |
|-------|
| |
| |

| BEB046100 (Single) | BEB046150 (King Single) | |
|---|----------------------------|--|
| 2288 mm | 2288 mm | |
| 900 mm | 1060 mm | |
| 110 - 740 mm | 110 - 740 mm | |
| 0 - 60° | 0 - 60° | |
| 0 - 35° | 0 - 35° | |
| Adjustable | Adjustable | |
| 250 kg | 250 kg | |
| 10 Years Frame, 3 Years Electronics, 1 Year Handset & Cables | | |

| Code | Description |
|-----------|--|
| BEB046200 | Tahoe Bed End Set - Sublime Teak - Single |
| BEB046220 | Tahoe Bed End Set - Sublime Teak - King Single |
| BEB046520 | Grab Assist Bar / Pole |
| BEB046525 | Split Folding Side Rails (Pair) |
| BEB046530 | IV Pole |
| BEB046540 | Self Help Pole |



Stock Bed Ends *BEB046200 (Tahoe) *BEB046220



Grab Assist Bar / Pole (BEB046520)



IV Pole (BEB046530)



Split Folding Side Rails (Pair) (BEB046525)



Self Help Pole (BEB046540)

Home Care & Hospital Beds



Malvestio Botero Bariatric Ward Bed - 400kg



Antibacterial frame coating for improved infection control











- 4 x section electrically operated bariatric ward bed
- Adjustable mattress platform with 102-122cm tool free adjustments
- Leg-rest adjustment is ergonomically synchronised with the knee break
- Integrated bed length extension up to 18cm
- IV & Self-Help Pole mounting brackets
- Electric height adjustment with 2 x heavy duty actuators
- Attendant control panel with intuitive button placement for bed control
- Provided battery backup as standard, allowing for power supply when transporting around wards or in the event of a power failure
- Translating backrest (Auto-Regression) shifts backwards by 11cm to compensate for patient movement towards the foot
- Folding side rails as standard
- Base frame is equipped with 150mm twin wheel castors featuring a centralised locking system that can be operated by foot pedals on either side of the bed end

| Code | BEB050200 | |
|------------------------|----------------|--|
| Overall Length | 2260 - 2440 mm | |
| Overall Width | 1020 - 1220 mm | |
| Overall Height | 380 - 850 mm | |
| Patient Surface Length | 2000 - 2180 mm | |
| Patient Surface Width | 1020 - 1220 mm | |
| Backrest Angle | 0° to 70° | |
| Knee Bend Angle | 0° to 30° | |
| Trendelenburg | -17° to 17° | |
| Castor Size | 150 mm | |
| SWL | 400kg | |
| Warranty | 2 Years | |

| Code | Description |
|-----------|---------------------------------|
| BEA825100 | Self Help Pole |
| BEA825101 | IV Pole |
| BEA825102 | Weigh Scale (Factory Fitted) |
| BEA825103 | Catheter Bag Holder |
| BEA825104 | Accessory Holder Bar |
| BEA825105 | Oxygen Bottle Holder |

Optional Accessories







Pole (BEA825101)



Malvestio Smile Paediatric Mattress - 150kg







Four layer foam, designed for larger users with a maximum user weight of 300kg. Full size memory foam overlay and side walls for added stability and ease of transfers. European fabric cover with fully welded seams, water resistant zippers and reinforced waterfall flaps for optimum infection control.

- Three function electrically adjustable hospital bed designed with growing people in mind
- Colourful side rails and decals to meet interior design themes and add to a more cheerful patient environment
- Capacity to lengthen at the press of a button to transition from child to young adult size proportions
- Biocote antibacterial surface treatment for improved infection control
- Durable blow molded side rails with gas spring dampeners for smooth operation
- Onboard battery provides ward mobility and emergency backup

- Integrated control panel in foot end section for carer and parent access
- Backrest anlge 0 to 70 degrees
- Knee bend angle 0 to 36 degrees
- Trendelenburg -15 to 15 degrees

| Code | PED490110 |
|----------------|-----------|
| Overall Length | 1580 mm |
| Overall Width | 790 mm |
| Overall Height | 140 mm |
| S.W.L | 150 kg |
| Warranty | 2 Years |

Malvestio Botero Static Mattress - 400kg







Four layer foam, designed for larger users with a maximum user weight of 300kg. Full size memory foam overlay and side walls for added stability and ease of transfers. European fabric cover with fully welded seams, water resistant zippers and reinforced waterfall flaps for optimum infection control.

- Two layer foam design with high density base and castellated surface caters to bariatric patients
- Multi-stretch cover with welded seams for stringent infection control
- Strong high density side walls improve patient transfer efficiency while cut back hinge points allow the mattress to easily move with the mattress platform

| Code | BEM050200 |
|----------------|-----------------------------|
| Overall Length | 2000 mm |
| Overall Width | 1220 mm |
| Overall Height | 200 mm |
| S.W.L | 400 kg |
| Warranty | 5 Years Foam, 2 Years Cover |

Mattress Overlays & Replacements













Aspire Lifecomfort Bariatric Foam Mattress - King Single



Four layer foam, designed for larger users with a maximum user weight of 300kg. Full size memory foam overlay and side walls for added stability and ease of transfers. European fabric cover with fully welded seams, water resistant zippers and reinforced waterfall flaps for optimum infection control.

- Pressure reducing bariatric mattress
- Extra weight capacity up to 300 kg
- Four layer foam support 200 mm high for bariatric use with maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress
- Two way waterproof stretch cover with welded seams
- PVC bottom cover
- Four sided zip for easy cover removal

- Stitched zipper cover with waterproof flap for excellent hygiene
- Independently tested and certified for durability and weight rating

| Code | BEM046700 (King Single) |
|----------------|-----------------------------|
| Overall Length | 2030 mm |
| Overall Width | 1050 mm |
| Overall Height | 200 mm |
| S.W.L | 300 kg |
| Warranty | 5 Years Foam, 2 Years Cover |










Aspire Lifecomfort Acute Care Foam Mattress

- Strong U-Shaped Foam Outer Structure Promotes Ease of Patient Transfers
- Welded High Temp Washable Fabric for Optimal Infection Control
- Vapour Permeable Four Way Stretch Fabric Improves Skin Microclimate
- Castellated Patient Surface for Improved Air Flow & Efficient Pressure Relief



| Code | BEM046615 (Single) | BEM046616 (King Single) |
|-----------------------|-----------------------------|----------------------------|
| Overall Length | 1980 mm | 2030 mm |
| Overall Width | 880 mm | 1050 mm |
| Overall Height | 150 mm | 150 mm |
| S.W.L | 250 kg | |
| Warranty | 5 Years Foam, 2 Years Cover | |
| Incl. Evac | BEM046615EV | BEM046616EV |









Aspire Lifecomfort Premier Care Acute Foam Mattress

- Pressure care mattress suitable for high risk patients
- Four layer foam support for maximum comfort and support
- Full size memory foam top layer with MEMair[™] Breeze open cell foam for improved microclimate
- Integrated cotton stockinette promotes foam longevity
- Reinforced side walls for patient stability and ease of transfers
- · Inbuilt handles for safe manual handling
- Two way waterproof stretch top cover
- Heavy duty PVC base cover



- Three sided zipper for easy cover
- Stitched zipper cover with waterproof flap for improved infection control
- PUPPS 2 (2004) Compliant
 - *Also available with Evacuation Sheet

| Code | BEM046735 (Single) | BEM046740 (King Single) |
|-----------------------|-----------------------------|----------------------------|
| Overall Length | 2000 mm | 2030 mm |
| Overall Width | 900 mm | 1050 mm |
| Overall Height | 150 mm | 150 mm |
| S.W.L | 250 kg | |
| Warranty | 5 Years Foam, 2 Years Cover | |
| Incl. Evac Sheet | BEM046735EV | BEM046740EV |

^{*}Also available with Evacuation Sheet





1300 133 120 www.aidacare.com.au



CA35.6 190719