

PRODUCT UPDATE



4H400E Electric Shower Trolley UPDATE!

The 4H400E Electric Shower Trolley supports efficient and safe showering and hygiene care of immobile patients. Electrically adjustable positioning promotes correct working posture, while side supports and soft padded liner provide comfort and support for the patient.



BTS147100 – 4H400E Electric Shower Trolley

Electric Adjustable Height: Supports ergonomic working positions for carers, through a 400mm range controlled by the hand control.

Stainless Steel Frame: Durable construction enhances product lifespan and reliability.

Soft Shower Liner: Provides patient comfort during showering and care.

Heavy Duty Castors: Three total-lock and one directional 125mm castors for controlled movement and stability.

Removable Battery: Allows fast battery change for continuous operation.