



BATH TRANSFER BENCH WITH ROTATING SEAT

The rotating bath transfer bench helps make bathing safer and more comfortable for people with reduced mobility. The seat is shaped to provide comfort during transfers, and its lightweight design makes it easy to move. The seat rotates 360°, allowing for smooth and secure transfers in and out of the bath, which reduces strain on both the user and caregiver.



- Comfortable padded backrest, seat and arms offer a warmer sitting surface
- Height adjustable legs with 12.5mm increments suits a range of bath to floor heights
- Contoured seat with a safety belt, slides smoothly and is easy to clean
- 360° rotating seat (with lockable positions every 90°) suits left and right bath orientations
- Slip resistant suction rubber tips for added stability
- Lightweight aluminium construction for ease of use and portability

BATH TRANSFER BENCH	Overall Height mm	Overall Length mm	Overall Width mm	Seat Width mm	Seat Depth mm	Seat Height mm	Total Weight kg
BTS150310	795 - 905	475	950	425	345	505 - 615	8.1