



# GENESIS HD ADJUSTABLE SEAT WALKER

The Aspire Genesis HD is extra strong to support users weighing up to 204kg, yet it is still lightweight and easy to lift and manoeuvre due to the strengthened aluminium elliptical tube frame. It is highly adaptable and comfortable with adjustable seat and handle heights, padded seat and generously wide padded backrest. Rugged solid 8" wheels are ideal for both indoor and outdoor terrain.

**Adjustable Seat Height:** Suitable for variable user sizes, optimising seat height for the user assists with sitting and standing easily and safely

**Adjustable Handle Height and Soft-Touch Brakes:** Correct handle height promotes walking posture while soft brakes are nice and easy to operate

**Wide Padded Seat & High Weight Capacity:**

Ensures comfort and stability for larger users, with the added benefit of a discreet storage pocket beneath the seat

**Large Padded Backrest:** Provides generous and comfortable lumbar support when seated

**Oversized Under-Seat Bag:** Ample storage space for personal items, accessible by lifting the seat

**Easy Folding Mechanism and Lightweight:** Lift to fold for compact storage and easy handling at only 11.3kg total weight

**Solid 8" Wheels with Inline Rear Brakes:** Tough treaded wheels for a smooth ride in variable terrain, while inline rear brakes maximise manoeuvring space for the user



MARBLE  
BLUE



GENESIS HD  
ADJUSTABLE  
SEATWALKER

Overall  
Height

Overall  
Length

Overall  
Width

Seat  
Width

Seat  
Height

Width  
Between  
Arms

Castor  
Size

Product  
Weight

SWL

SWL  
Under  
Seat Bag

mm

mm

mm

mm

mm

mm

mm

kg

kg

kg

WAF750336

875 -  
1050

700 - 910

660

450

500-590

520

200

11.3

204

5

\*Range of motion dimensions are accurate to a +/-10mm tolerance