

# HS118 ULTRALIGHT COMPACT SCOOTER

Aluminium Frame Pulls Apart Easily, Ready for Any Journey!

The HS118 is an ultralight portable scooter, designed for easy lifting in and out of a vehicle. Packed with category leading features, this scooter offers a 14km range and is perfect for indoor and outdoor use.

Heaviest Part  
ONLY  
**10.2kg**



- Lightweight aluminium frame extends travel range
- Disassembles into 5 pieces with the heaviest component weighing only 10.2kg
- Max speed of 6.4km/h with a max. range of 14km
- Adjustable steering to achieve your perfect driving position
- Swivelling seat with adjustable armrests make entry, exit and transfer easy!

**RED**

**2 YEAR WARRANTY**  
Frame, Motors & Brakes

**1 YEAR WARRANTY**  
Battery

**SWL 115KG**

SCT671115	Overall Length	Overall Width	Overall Height	Seat Width	Wheel Size	Max. Speed	Ground Clearance	Turning Radius	Travel Range	Battery Size	Product Weight incl. Batteries
	mm	mm	mm	mm	mm/in	km/h	mm	mm	km	V	kg
	1010	550	890	430	195/8	6.4	35	1350	14	1 x 11.6AH	31