



# MAUI DUAL ACTION 2 MOTOR

- A luxurious lift recline Chair, the Maui Dual Action delivers multiple positioning options, with additional backrest padding for a plush sitting and reclined experience
- Finished with soft quilted upholstery
- Dual motor lift recline chair
- Memory foam seating surface – improves comfort and helps to facilitate pressure redistribution
- Independent control of back and legrest functions
- Head and armrest covers for improved longevity
- Available in mocha quilted upholstery



Back Recline



Leg Recline



Colour



Mocha  
Quilted Upholstery

Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Mocha – CHP227110	500mm	500mm	540mm	630mm	1100mm	780mm	130kg

ASBR26 210623