

SIGNATURE DUAL ACTION

2 MOTOR



Mink Warwick
Global Microsuede



Lagoon Warwick
Global Microsuede



Heritage Vinyl



Available in
custom upholstery



Optional Neck Pillow

CHP197781 | Mink Warwick
CHP197783 | Lagoon Warwick
CHP197782 | Heritage



Back Recline



Leg Recline

- Available in three sizes of single and dual motor models
- Dual action recliner – dual control of leg and recline for infinite positioning
- Depth adjustable seat – a choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning
- Easy to move – four swivel castors facilitate movement. Provided as standard (seat height increases when fitted)
- Features an all-day battery backup to maximise safety and independence
- Comes with integrated lithium battery backup

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228200	Mink Warwick	Small	480mm	460 – 385mm	480mm	580mm	980mm	740mm	155kg
CHP228210	Lagoon Warwick								
CHP228220	Heritage								
CHP228230	Mink Warwick	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	155kg
CHP228240	Lagoon Warwick								
CHP228250	Heritage								
CHP228260	Mink Warwick	Large	500mm	510 – 535mm	540mm	650mm	1110mm	780mm	155kg
CHP228270	Lagoon Warwick								
CHP228280	Heritage								