

MAXI SHOWER CHAIR











- Height adjustable legs with easy C-Clip adjustment points for changing needs
- Heavy duty reinforced frame that accommodates very high SWL
- Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
- · All frames individually treated for corrosion resistance
- Contoured plastic seat with drainage holes and windows to provide postural support, improve drainage and reduce skin adhesion
- Additional height increments on one leg to accommodate uneven surfaces
- Raised armrests for assistance in transferring on and off the stool







MAXI shower chair	Overall Height	Overall Width	Overall Depth	Seat Height	Seat Width	Seat Depth	Between Arms
	mm	mm	mm	mm	mm	mm	mm
BTS118020	670 - 1005	645 - 845	495	455 - 575	575 - 775	455	575 - 775



DVA CONTRACTED



AS / NZS Approved



MAXI is registered in the