



# IDAHO DUAL ACTION 2 MOTOR

- Contemporary medium size
- Gentle tapered armrests suit the modern home environment
- Simultaneous adjustable back and leg support functions
- Integrated lift and recline functions for easy sit-to-stand
- Generous padding in backrest for improved lumbar support



Leg Recline



Back Recline



**Optional Neck Pillow**

CHP227171 | Charcoal Fabric

CHP227172 | Latte Vinyl



Charcoal Fabric



Latte Vinyl

Code	Colour	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP227160	Charcoal	480mm	480mm	480mm	620mm	1100mm	800mm	150kg
CHP227165	Latte							

# SIGNATURE SINGLE ACTION 1 MOTOR SPACE SAVER



### Optional Neck Pillow

CHP197781 | Mink Warwick  
CHP197783 | Lagoon Warwick  
CHP197782 | Heritage



Full Recline Space Saver



Mink Warwick  
Global Microsuede



Lagoon Warwick  
Global Microsuede



Heritage Vinyl



Available in  
custom upholstery

- Available in three sizes of single and dual motor models
- Space saver recliner – space saving recline for close wall positioning
- Single motor – back and legs move simultaneously
- Features and all-day battery backup to maximise safety and independence
- Depth adjustable seat – a choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning
- Four swivel castors facilitate movement. Provided as standard (seat height increases when fitted)
- Comes with integrated lithium battery backup

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228100	Mink Warwick	Small	480mm	460 – 485mm	480mm	580mm	980mm	740mm	145kg
CHP228110	Lagoon Warwick								
CHP228120	Heritage								
CHP228130	Mink Warwick	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	145kg
CHP228140	Lagoon Warwick								
CHP228150	Heritage								
CHP228160	Mink Warwick	Large	500mm	510 – 535mm	540mm	650mm	1110mm	780mm	145kg
CHP228170	Lagoon Warwick								
CHP228180	Heritage								

# SIGNATURE DUAL ACTION

## 2 MOTOR



Mink Warwick  
Global Microsuede



Lagoon Warwick  
Global Microsuede



Heritage Vinyl



Available in  
custom upholstery



**Optional Neck Pillow**

CHP197781 | Mink Warwick  
CHP197783 | Lagoon Warwick  
CHP197782 | Heritage



**Back Recline**



**Leg Recline**

- Available in three sizes of single and dual motor models
- Dual action recliner – dual control of leg and recline for infinite positioning
- Depth adjustable seat – a choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning
- Easy to move – four swivel castors facilitate movement. Provided as standard (seat height increases when fitted)
- Features an all-day battery backup to maximise safety and independence
- Comes with integrated lithium battery backup

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228200	Mink Warwick	Small	480mm	460 – 385mm	480mm	580mm	980mm	740mm	155kg
CHP228210	Lagoon Warwick								
CHP228220	Heritage								
CHP228230	Mink Warwick	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	155kg
CHP228240	Lagoon Warwick								
CHP228250	Heritage								
CHP228260	Mink Warwick	Large	500mm	510 – 535mm	540mm	650mm	1110mm	780mm	155kg
CHP228270	Lagoon Warwick								
CHP228280	Heritage								