



# NINJA CARBON FIBRE SEAT WALKER

The Aspire Ninja is a lightweight seat walker designed for everyday use. Made from carbon fibre, it's easy to lift in and out of a car and has large front wheels for greater mobility outdoors.



Minimum Weight  
**ONLY 5.4kg**



### Compact Wheelbase

Rear wheels are positioned for maximum clearance during entry and exit from the frame



### Easy Folding

The walker folds easily for convenient transport and compact storage



**5** YEARS WARRANTY  
*Frame*

**1** YEAR WARRANTY  
*Components*



**Great Outdoors**  
The large 10" wheels provide stability and smooth movement across a variety of terrains



**Adjustable Handles**  
Allows for personalised positioning, helping to maintain proper posture



**Soft-Touch Brakes**  
Lock the rear wheels to enhance safety and stability during use



**Utility Bag**  
Holds up to 5kg, offering convenient storage for personal items



**Ergonomic Handles**  
These offer a comfortable grip during extended use. The flat-top design provides a stable resting place for arms while seated



**Compact Design**  
Improves turning, offering smoother navigation in tight spaces

## Discover Our Exclusive Range of Accessories



**Adjustable Height IV Pole**  
WAA695650



**Oxygen Bottle Holder Bag**  
WAA695600



**Night Light LED Lamp**  
WAA695630



**Crutch & Cane Holder - comes as standard**

## NINJA CARBON WALKER

	Overall Height	Overall Length	Overall Width	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**
	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	mm/in	kg/st	kg/st
WAF705380 MEDIUM	860 - 945 / 34 - 37	710 / 28	625 / 24.5	460 / 18	540 / 21	245 / 9.6	F - 250 / 10 R - 200 / 8	5.4 / 12	6.35 / 14
WAF705390 TALL	860 - 1025 / 34 - 37	710 / 28	625 / 24.5	460 / 18	620 / 24.4	245 / 9.6	F - 250 / 10 R - 200 / 8	5.5 / 12	6.4 / 14

\*Range of motion dimensions are accurate to a +/-10mm tolerance  
\*\*Includes accessories, bag, and back strap.