

# aspire PACE WALKING TUTOR

The Aspire Pace Walking Tutor is designed to support balance and stability while helping the user to walk with an upright gait. This unique aid is ideal for individuals with limited hand and wrist strength as lower arm supports allow them to bear upper body weight comfortably and confidently.

## FEATURES & BENEFITS



### Comfortable Handles and Lower Arm Supports

Ergonomic Handles with Brakes  
Depth and angle adjustable to ensure optimum comfort and accessibility to brakes and handles for control and support.



### Height-Adjustable Design

The soft padded lower arm support platform (including brakes & handles) has lever activated height adjustment to achieve the optimum height for the individual at any time, and to adjust as needs change. Includes gas assisted raise and lockable height setting once adjusted.



### Efficient Braking System

Ensures reliable performance while providing intuitive rear-wheel control, requiring minimal maintenance.



### Swivel Castors

All castors are total locking to ensure the tutor can remain safely stationary and stable when required



**1** YEAR WARRANTY



Luxury Look & Feel  
Danish Designed



Overall Length



Overall Height



Overall Width



Width Between Forearm Supports



Wheel Size



Product Weight



Pace Walking Tutor  
WAF699220

800mm

915 – 1325mm

640mm

660 – 790mm

100mm

18kg

31.5"

36 – 52"

25"

26 – 31"

4"

39.7lbs

Australia  
1300 133 120 | [aspirecare.com.au](http://aspirecare.com.au)

United Kingdom  
01622 541235 | [aidacare.co.uk](http://aidacare.co.uk)