




# SHOWER CHAIR WITH ARMS



- Height adjustable legs with easy C-Clip adjustment points for changing needs
- Reinforced frame that accommodates a high SWL
- Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
- All frames individually treated for corrosion resistance

- Contoured plastic seat with drainage holes and windows to provide postural support, improve drainage and reduce skin adhesion
- Additional height increments on one leg to accommodate uneven surfaces
- Raised armrests for assistance in transferring on and off the stool

**5** YEARS WARRANTY  
Corrosion

**1** YEAR WARRANTY  
Structure



## SHOWER CHAIR WITH ARMS

	Overall Height mm	Overall Width mm	Overall Depth mm	Seat Height mm	Seat Width mm	Seat Depth mm	Between Arms mm	Product Weight kg
STEEL BTS118000	735 - 910	512	570	375 - 550	440	400	470	3.4
ALUMINIUM BTS118005	735 - 910	512	570	375 - 550	440	400	470	3.1
WIDE - STEEL BTS118010	735 - 910	570	570	375 - 550	440	400	525	3.9
WIDE - ALUMINIUM BTS118015	735 - 910	570	570	375 - 550	440	400	525	3.3

