aspire

TOILET SURROUND ALUMINIUM

FEATURES & BENEFITS



Raised armrests for assistance transferring on and off a toilet



Height adjustable legs with easy push button height adjustments. Additional increments also available on one leg to accommodate uneven surfaces



Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability



TOILET SURROUND	Overall Height	Overall Width	Overall Depth	Between Arms	Product Weight
	mm	mm	mm	mm	kg
BTT155100	535 - 660	500	435	450	3



DVA CONTRACTED



AS / NZS Approved

