



VOGUE FOREARM WALKER


For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:



White / Black



1 YEAR WARRANTY



Padded Forearm Support



Folds to 8.8"

✓ AS / NZS Approved

✓ Vogue is registered in the ARTG, AUST R ARTG 271255

FEATURES & BENEFITS

- **Ergonomic Design** offers personalised height adjustment which promotes correct posture
- **Handbrakes** lock for user safety and stability
- **Simple** folding mechanism locks for easy transportation and storage
- **Storage** removable bag carries up to 5kg of essentials
- **Bonus** walking cane holder



FOREARM WALKER	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
	mm /in	mm /in	mm /in	mm /in	mm /in	mm /in	mm /in	mm /in	kg /lbs	kg /lbs	kg /lbs
WHITE/BLACK	655	608	1020-1270	240	441	620	225	203	9.1	10	5
WAF705400	26	24	40-50	9.5	17.5	24.5	9	8	20	22	11

** Includes accessories, bag and backstrap.

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk