

LIGHTWEIGHT 2 SEAT WALKER

This modern European designed seat walker has a clean and crisp frame design and is available in a selection of eye-catching colours and prints to really stand-out while keeping you on the move!

**UNDER
6KG**



FEATURES AND BENEFITS

- Lightweight aluminium body improves manoeuvrability and portability
- Ergonomic handles with personalised height adjustment to promote comfortable and correct posture
- Compact turning circle for improved mobility in tight spaces
- Simple folding mechanism locks for easy transportation and storage
- Handbrakes lock for user safety and stability
- Removable bag carries up to 5 kg of essentials
- Deep seat and wide backrest for improved comfort and support
- Bonus walking cane holder



ZEBRA

WAF705310ZB

1 YEAR WARRANTY

SWL
150kg



LEOPARD

WAF705310LE



FLOWER

WAF705310FL



RED

WAF705310RE



CHAMPAGNE

WAF705310CH



Overall Length
mm

655



Overall Width
mm

620



Overall Height
mm

785 – 920



Seat Depth
mm

240



Seat Width
mm

465



Seat Height
mm

540



Folded Width
mm

220



Wheel Size
in

8"



Minimum Functional Weight
kg

5.4



Typical Functional Weight**
kg

6.45



Bag SWL
kg

5

** Includes accessories, bag and backstrap.