



WHEELCHAIR ACCESSORIES

Aspire for... Comfort, stability, support



aspire WHEELCHAIR ACCESSORIES



Anti-tip Bars

Ensures chair does not tilt over when traveling across inclines and uneven surfaces.

AVAILABLE FOR

Vida | Vida X | Vida MG | Sociaite | Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward | Rehab RX 2 | Rehab RX 2 JNR



Attendant Brake Kit

Kit includes 12" wheels and attendant brakes, or attendant brakes for 24" self-propelled wheels.

AVAILABLE FOR

Evoke 2



Safety Lap Belt

Suitable for supporting users with low tone or weakness. Metal push button for secure fit.

AVAILABLE FOR

Vida | Vida X | Vida MG | Sociaite | Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward | Rehab RX 2 | Rehab RX 2 JNR



Contoured Padded Overlay

Offers increased comfort, ideal for users in their wheelchair for longer periods of time. Breathable washable cover.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward



Drum Brake Kit

This upgrade kit can be retrofitted, providing the attendant with more control on sloped surfaces.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Rehab RX 2 | Rehab RX 2 JNR



Oxygen Bottle Holder

Holds oxygen bottle in place. Designed to hold C size bottles.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Rehab RX 2 | Rehab RX 2 JNR | Transit 3 | Assist 3



Stump Support

For amputee support.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Rehab RX 2 | Rehab RX 2 JNR | Transit 3 | Assist 3



Elevating Legrest

Padded legrests for improved comfort.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward | Rehab RX 2 | Rehab RX 2 JNR



IV Pole

Facilitates intravenous infusion.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward | Rehab RX 2 | Rehab RX 2 JNR



Extension Lever

Extends brake lever for user accessibility.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Assist 3 Ward



Calf Strap

Provides extra support for user's calves.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Transit 2 | Assist 2 | Transit 3 | Assist 3 | Asisst 3 Ward | Rehab RX 2 | Rehab RX 2 JNR



Base Frame Extension Kits

Extends stump support to cater for taller users.

AVAILABLE FOR

Rehab RX 2 | Rehab RX 2 JNR



Depth & Height Adjustable Armrests

Provides multiple support positions for users arms.

AVAILABLE FOR

Evoke 2



Extended Back Canes

Extends back canes up to 50mm to facilitate taller carers.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR



Stroller Handle

Angle adjustable handle allows greater ease of use for a carer.

AVAILABLE FOR

Rehab RX 2 | Rehab RX 2 JNR



Lateral Chest Support

Facilitates postural positioning.

AVAILABLE FOR

Evoke 2 | Evoke 2 JNR | Rehab RX 2



Deluxe Headrest

Angle adjustable head support. Helps to promote user posture and comfort.

AVAILABLE FOR

Evoke 2




Gutter Armrest Pad

Provides extra support for users with weak or unsupported arms.

AVAILABLE FOR

Evoke 2



Adjustment Kit

Allows for easy, tool-free adjustment of a wheelchair.

AVAILABLE FOR

Rehab RX 2 | Rehab RX 2 JNR



Meal Tray

Provides sturdy surface for meal time.

AVAILABLE FOR

Rehab RX 2

Aspire Rehab RX2 Wheelchair Accessories

○ Standard Feature Included ● Compatible

Code	Description	Rehab RX2	Rehab RX2 JNR
Stump Supports			
WCA848220	Left - 300/320/350/400/410mm	●	●
WCA848230	Right - 300/320/350/400/410mm	●	●
WCA848240	Left - 500/550/570/610mm	●	
WCA848250	Right - 500/550/570/610mm	●	
Lateral Chest Support Assembly			
WCA848450	Left	●	
WCA848460	Right	●	
Elevating Legrests			
WCA848100	Left - 300/320/350/400/410mm	●	●
WCA848110	Right - 300/320/350/400/410mm	●	●
WCA848120	Left - 450/460/500/510mm	●	
WCA848130	Right - 450/460/500/510mm	●	
Oxygen Bottle Holders			
WAA695600	Oxygen Bottle Bag	●	●

Code	Description	Rehab RX2	Rehab RX2 JNR
Rear Wheels with Drum Brake Kits			
WCA848510	16"	●	●
WCA848520	20"	●	○
WCA848500	24"	○	○
Base Frame Extension Kits			
WCA848530	50mm - Black	●	
WCA848540	100mm - Black	●	
Other			
WCA848550	Tool-free Adjustment Kit	●	●
WCA848570	Meal Tray - 460mm to 510mm	●	
WCA848270	Anti-tip Bars	○	○
WCA848280	IV Pole	●	●
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	●	●
WCA848580	Calf Strap	●	●
WCA848590	Stroller Handle	●	●

*Note Side transfers not possible when gutter armrests are attached

Aspire Vida Range Wheelchair Accessories

Code	Description
WCA848275	Anti-tip Bars
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle

Aspire Socialite Wheelchair Accessories

Code	Description
WCA848270	Anti-tip Bars
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle

Aspire Evoke 2, Transit 2, Assist 2, Transit 3 & Assist 3 Wheelchair Accessories

Code	Description	Evoke 2	Evoke 2 JNR	Transit 2	Assist 2	Transit 3	Assist 3
Contoured Padded Overlays							
WCA848325	300mm		●				
WCA848330	350mm	●					
WCA848340	400mm	●		●	●	●	●
WCA848350	450mm	●		●	●	●	●
WCA848360	500mm	●			●	●	●
WCA848370	570mm	●			●	●	●
WCA848380	610mm	●			●	●	●
Stump Supports							
WCA848220	Left - 300 to 410mm	●	●				
WCA848230	Right - 300 to 410mm	●	●				
WCA848240	Left - 500 to 610mm	●					
WCA848250	Right - 500 to 610mm	●					
WCA848260	Extension Bracket	●	●				
Lateral Chest Support Assembly							
WCA848430	Left	●	●				
WCA848440	Right	●	●				
Elevating Legrests							
WCA848100	Left - 300 to 410mm	●	●	●	●		
WCA500080	Left - 400mm					●	●
WCA848110	Right - 300 to 410mm	●	●	●	●		
WCA500085	Right - 400mm					●	●
WCA848120	Left - 450 to 510mm	●		●	●		
WCA500070	Left - 450 to 610mm					●	●
WCA848130	Right - 450 to 510mm	●		●	●		
WCA500075	Right - 450 to 610mm					●	●
WCA848140	Left - 570/610mm	●			●		
WCA848150	Right - 570/610mm	●			●		
WCA848080	Left - Evoke 2 JNR		●				
WCA848090	Right - Evoke 2 JNR		●				
Rear Wheels with Drum Brake Kits							
WCA848510	16"		●				
WCA848520	20"		●				
WCA848500	24"	●					

Code	Description	Evoke 2	Evoke 2 JNR	Transit 2	Assist 2	Transit 3	Assist 3
Depth & Height Adjustable Armrests							
WCA848390	Left	●					
WCA848400	Right	●					
Gutter Armrest Pads*							
WCA848410	Left	●					
WCA848420	Right	●					
Brake Extension Levers							
WCA848310	Left	●	●	●	●	●	●
WCA848320	Right	●	●	●	●	●	●
Brake Kits							
WCA848592	Attendant Brakes	●					
WCA848594	12" Wheel & Brake Kit	●					
Other							
WCA848290	Oxygen Bottle Holder	●	●				
WCA848270	Anti-tip Bars with Mounting Hardware - pair	●	●	●	●		
WCA500050	Anti-tip Bars with Mounting Hardware - pair					●	●
WCA848272	Advanced Depth-adjustable Anti-tip Bar	●					
WCA848280	IV Pole	●	●	●	●		
WCA848285	IV Pole with Mounting Clamp					●	●
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	●	●	●	●		
WCA500055	Safety Lap Belt - Metal Seatbelt Buckle					●	●
WCA848470	Extended Back Cane - Pair - 19"/21"	●	●				
WCA848480	Deluxe Headrest Assembly 400 to 450mm	●	●				
WCA848481	Deluxe Headrest Assembly 500 to 610mm	●					
WCA848580	Calf Strap	●	●	●	●		
WCA500065	Calf Strap					●	●
WCA500060	Tool Kit					●	●
AAL848590	Security Mod			●	●	●	●

*Side transfers not possible when gutter armrests are attached

