

# ALTITUDE CHAIR

Aspire for... Outcomes, Comfort and Mobility



The **ALTITUDE Vertical Lift** chair is designed to accommodate varying user and facility requirements for enhanced comfort and efficient care provision.



▲ Assisted transfer functionality

## FEATURES & BENEFITS

- **Mobile patient care:** Integrated power source allows full chair functionality throughout the care environment for all-day independence.
- **Comfort and pressure relief:** Multi-stretch vapour permeable seat and backrest help to reduce shear and promote optimal skin microclimate. The outer chair surface is upholstered in medical-grade vinyl for longevity and infection control.
- **Intuitive control:** Large backlit buttons make operation a simple process and enhance user independence.
- **Assisted transfer functionality:** Vertical lift and removable armrests facilitates safe sit-to-stand and lateral transfers, simplifying patient handling.
- **Oncology armrests (CHP227710):** Designed to enhance patient comfort, the wide forearm support assists with effective administration of treatment such as blood transfusions.
- **Available in Bitter Chocolate (multi-stretch fabric) and Cream (vinyl).**

**1 YEAR WARRANTY**  
BATTERY + HAND CONTROL

**2 YEARS WARRANTY**  
FABRIC + MOTOR + MECHANISM

**S.W.L 180KG**

CHP227700 Aspire ALTITUDE Vertical Lift chair	SEAT HEIGHT	SEAT DEPTH	SEAT WIDTH	SEAT TO TOP OF BACK	OVERALL HEIGHT	OVERALL DEPTH	OVERALL WIDTH	OVERALL RECLINE
	560 - 700mm	520mm	510mm	710mm	1200 - 1340mm	800mm	710mm	1780mm

ASFY617\_300120

1300 133 120

www.aspirecare.com.au

