ASPIRE ALUMINIUM COMMODE

Mobility | Support | Hygiene





Soft Opening PU Seat



Swing Away Footrests



24" Lightweight Wheels (Self Propelled model)

FEATURES & BENEFITS

- Ideal for in home and community use
- Lightweight Aluminium design
- Easy clean, padded seat, back and arms
- Tool free height adjustment allows the commode to be easily adapted to varying needs
- Swing back comfort arms assist transfers
- Bowl and carrier for improved infection control
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- 4 x Locking castors on the attendant propelled model





	OVERALL WIDTH (mm)	OVERALL HEIGHT (mm)	SEAT WIDTH (mm)	BACKREST HEIGHT (mm)	ARMREST HEIGHT (mm)	ANTERIOR SEAT HEIGHT (mm)	POSTERIOR SEAT HEIGHT (mm)
Aspire Aluminium Commode				1	Ī		
Attendant Propelled BTC068010	460	989 - 1039	440	495	220	575	580
Self Propelled BTC068020	610	969 - 994	440	490	235	575	580





