

# NEW PRODUCT LAUNCH



## Aspire Lifestyle Rotating Bed – New Product RELEASE!

This bed supports diverse care needs with advanced profiling features, including rotation and tilt functions. It smoothly transitions between lying, sitting, and standing, improving comfort for both users and caregivers. Ideal for home, nursing, and residential care settings.



HI LO



BACK RAISE



LEG RAISE



SEATED



TILT TO EXIT

- Aspire Lifestyle Rotating Bed – Single – BE046045 – SWL 220kg
- Aspire Lifestyle Rotating Bed – Mattress – BEM046300 – SWL 120kg
- Aspire Lifestyle Rotating Bed – HD Mattress – BEM046302 – SWL 180kg

# NEW PRODUCT LAUNCH

## PRODUCT FEATURES

- **Platform Rotation:** The upper section of the bed supports users during transitions from lying to sitting to standing. Its chair-like motion makes transfers smoother and allows caregivers to position themselves closer for safer assistance. The outer edges of the backrest platform wrap around the user before standing for additional upright support against gravity.
- **Low Height Adjustment:** Adjustable from 345mm, the bed supports sit-to-stand transitions for shorter users and can be raised with the press of a button to accommodate taller individuals.
- **Timber End Panels:** Stylish timber panels enhance the appearance, making it ideal for homecare settings.
- **Transportable Design:** Includes tool-free assembly and standard mounting trolleys, enabling quick preparation for delivery, transportation, and storage.
- **Central Locking Castors:** Easily lock all wheels with the touch of a foot, immobilising the bed in preparation for safe transfers.
- **Durable Side Rails:** Timber laminate rails adjust to multiple heights and lower as needed. These comply with IEC 60601-2-52 standards for safe and reliable care.
- **Tool-Free Setup:** Quick installation without removing electronics, ensuring efficient and hassle-free assembly.