

## **AFFINITY GRAB RAIL – Now Available!**

The Affinity Grab Handle/Bed Rail provides reliable and stable support for users during sitting, standing, and transitioning between lying and sitting positions.

The support frame fastens directly to the bed base to ensure the rail cannot shift during transfers which reduces falls risk.

The Affinity Grab Rail will now be available at all Aidacare locations. This successful product has come to us via David and Jackie at Mack'n'Me, former network partners, and now part of Aidacare!



Product Code: BEA013650



## NEW PRODUCT LAUNCH BENEFITS

- Support frames slide beneath the mattress and attach securely to the bed base, providing additional stability during use.
- Safety straps fasten beneath the bed to prevent shifting and eliminate trip hazards.
- Compatible with single and king single bed sizes.
- Rail can be installed on either side of the bed.
- Easily adjustable in height to accommodate a range of mattress thicknesses, ensuring a secure fit
- Includes a convenient accessory pouch for storing personal items.
- Provides sturdy and reliable support, making it easier to get in and out of bed.

