

NEW PRODUCT LAUNCH



AEROLITE CARBON FIBRE POWERCHAIR - NOW AVAILABLE!

The AeroLite Powerchair is light... VERY light. The modern carbon fibre frame is so light that it is a cinch to fold, lift and carry, to easily load into a car or to store when not in use.

AeroLite's lithium-ion batteries are approved for air travel and yes it comes with two of them! Perfect for taking you further for longer outings with no airport stress and less hassle and worry about battery charging downtime.

Intuitive joystick control, contoured comfort seating and puncture-proof tyres make for confident travel on different surfaces, indoors and out.

The AeroLite makes everyday life easier for users and carers alike.



Aspire AeroLite Carbon Fibre Powerchair

PWS654150

NEW PRODUCT LAUNCH

PRODUCT FEATURES

- **Side Folding Frame:** Stands up when folded, for ease of lifting into the car
- **Extended Range:** 2 x lithium-ion batteries provide a leading 20km range
- **Puncture Proof Tyres:** Provide confidence and reliability outdoors
- **Flip-up Footplates:** Improves transfer safety during seat entry and exit
- **Foldable Design:** Ideal for travel and easy to lift into the car
- **Tool-Free Assembly/Disassembly:** Quick and hassle-free setup
- **Ergonomic Seat Contouring:** Enhances comfort during prolonged sitting
- **Anti-Tippers:** Included as standard