

# NEW PRODUCT LAUNCH



## ASPIRE VOGUE NINJA WALKER AVAILABLE NOW!

Introducing the newest addition to our popular Carbon Fibre rollator family – the Aspire Vogue Ninja!

This new design features a compact wheelbase, maximizing clearance around the rear wheels for easier entry and exit. Larger front castors enhance performance on unsealed surfaces, while ergonomically shaped handles offer a comfortable arm position when resting.



WAF705380 - Aspire Ninja Seat Walker – Medium

WAF705390 - Aspire Ninja Seat Walker – Large

## NEW PRODUCT LAUNCH



### PRODUCT FEATURES

- **Compact Wheelbase:** Rear wheels are positioned for maximum clearance during entry and exit from the frame.
- **Compact Design:** Improves turning, offering smoother navigation in tight spaces.
- **Great for Outdoors:** The large 10" front wheels provide stability and smooth movement across various terrains.
- **Multiple Sizes:** Available in both Medium and Tall sizes.
- **Soft-Touch Brakes:** Lock the rear wheels to enhance safety and stability during use.
- **Ergonomic Handles:** Offers a comfortable grip during extended use. The flat-top design provides a stable resting place for arms while seated.
- **Height Adjustable Handles:** Allows for personalised positioning and correct posture.
- **Easy Folding:** Folds easily for convenient transport and compact storage.