

# NEW PRODUCT LAUNCH



## Aspire Vida MG Wheelchair RELEASE!

The VIDA MG, made from ultra-light magnesium and weighing just 6.7kg, opens up a world of freedom and convenience—perfect for those who love to travel and explore without limits.

The ergonomic contoured seat shape maximises comfort and supports posture by preventing users from sliding forward in the seat, perfect for a comfortable ride outdoors!



Minimum  
Chair Weight  
**6.7kg**

MWS449872 - Aspire VIDA MG Wheelchair - AP - 450mm - 113kg

MWS449873 - Aspire VIDA MG Wheelchair - SP - 450mm - 113kg

# NEW PRODUCT LAUNCH

## PRODUCT FEATURES



- Magnesium Frame - lightest in it's class, yet strong and durable.
- Intuitive folding frame and hinged backrest make transportation compact and easy!
- Ergonomic contoured seat shape provides superior comfort and postural support.
- Attendant operated hand brakes with parking lock
- Removable and adjustable leg rests for compact storage
- Rear storage pocket for essentials
- SWL of 113kg
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)

## DISTRIBUTION INFORMATION

- Stock to be held at all Distribution Centres and Showrooms.
- Stock has arrived in Moorebank and will be distributed to branches as per the planning setting in the coming week. Managers, please contact *Brett Heffernan* or *Duc Nguyen* to amend the planning settings (if needed).

## MARKETING AND PRICING

- Flyer available. ([link](#))
- Available on our website. ([link](#))
- List, Network Partner, Wholesale and External Sales Guide pricing available in GP.