


# PRODUCT UPDATE



## Aspire LIFECOMFORT Acute Care Mattress – Now with Heel Taper

Our Aspire LIFECOMFORT Acute Care Mattresses have undergone a design improvement and now include a gentle heel taper at no additional cost.

This improvement provides similar heel protection to the Aspire LIFECOMFORT Hybrid Comfort Mattress, redistributing pressure away from the vulnerable heel prominence to the larger surface area of the calf.



### Key Features

- Ultra-soft Heel Protection Zone improves pressure relief for vulnerable points. By gently tapering the heel area, pressure is redistributed effectively onto the less vulnerable calf area.
- Castellated Surface collapses in a wave-like movement reducing shear.
- Refer to the LIFECOMFORT Mattress brochure or website for full details and features.