

## Aspire Matisse Lift Chair in Petite Size UPDATE!

Due to internal feedback and market demand, we are pleased to announce the addition of a PETITE sized Aspire Matisse Quattro Lift Chair.

It retains all the same functionality and comfort features of the large Matisse, but with reduced seat depth, width and height ideal to support petite to small sized users.

The SWL is also reduced to 136kg to suit





## Sizing

Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Petite - CHP227555	455mm	455mm	510mm	560mm	1015mm	915mm	136kg
Large - CHP227560	500mm	545mm	570mm	730mm	1065mm	965mm	159kg



Issued: 22<sup>nd</sup> July 2024