

# PRODUCT UPDATE



## New Padded Seat for Aspire Over Toilet Aids

The padded toilet seat is an alternative to the standard white plastic toilet seats supplied with Aspire OTA's. It is 425mm wide, has a closed front and comes ready to clip-on.

### PRODUCT FEATURES

- **Padded Cushion:** Provides comfort for patients, particularly those with a low BMI or those needing extended sitting time.
- **Hygienic Design:** Easy to wipe clean, ensuring effective infection control.
- **Simple Installation:** Attaches to the frame with two small clips.

