

Litestream & iGo Power Wheelchair CHANGES!

The Pride Litestream Folding Carbon Fibre Power Wheelchair and the Pride iGo Folding Power Wheelchair now feature a redesigned joystick controller. While the controller's appearance has changed, its functionality remains the same.



- Power button: now positioned at the top.
- Speed increase button: located on the right.
- Speed decrease button: situated on the left.
- Horn button: placed at the bottom.
- Additional buttons: two at the bottom, currently without any function.

This change ensures a more intuitive layout, making it easier for users to access essential controls while maintaining familiar operation for existing users.

PRODUCT CODES

PWS654175 - iGo Folding Power Wheelchair

PWS654185 - Litestream Folding Carbon Fibre Electric Wheelchair



Issued: 15th November 2024