

NEW PRODUCT LAUNCH



Aspire Pace Walking Tutor – Now Available!

We are pleased to announce the release of the Aspire Pace Walking Tutor, our latest walking support solution that Replaces the discontinued Freedom Viking and Titan walkers.

The Aspire Pace Walking Tutor is designed to support balance and stability while helping the user to walk with an upright gait. This unique aid is ideal for individuals with limited hand and wrist strength as lower arm supports allow them to bear upper body weight comfortably and confidently.



FEATURES & BENEFITS

Ergonomic Handles with Brakes

Depth and angle adjustable to ensure optimum comfort and accessibility to brakes and handles for control and support.

Height-Adjustable Design

The soft padded lower arm support platform (including brakes & handles) has lever activated height adjustment to achieve the optimum height for the individual at any time, and to adjust as needs change. Includes gas assisted raise and lockable height setting once adjusted.

Intuitive Braking System

Easily activated brakes provide intuitive control of rear castors to assist support and balance while walking.

Swivel Castors

All castors are total locking to ensure the tutor can remain safely stationary and stable when required.

Safe Working Load: 150kg.