



4H400E ELECTRIC SHOWER TROLLEY

The 4H400E Electric Shower Trolley supports safe bathing and showering of immobile patients. It has features that improve care efficiency including adjustable sides, electrically adjustable positioning, and heavy-duty castors.



Electric Adjustable Height: 400mm range, controlled by the hand control supports ergonomic working positions for carers.

Stainless Steel Frame: Durable construction enhances product lifespan and reliability.

Soft Shower Liner: Provides patient comfort during showering and care.

Heavy Duty Castors: Three total-lock and one directional 125mm castors for controlled movement and stability.

Ergonomic Push Handle: Reduces strain and improves manoeuvrability for carers.

Removable Battery: Allows fast battery change for continuous operation.



1 YEAR WARRANTY
Components, Hand Control & Cables

5 YEAR WARRANTY
Frame

3 YEAR WARRANTY
Electronics



Max Total Load
(Patient & Water)

Code	Description	Total Length	Total Width	Total Height
BTS147100	4H400E Electric Shower Trolley	1890mm	730mm	780 - 1180mm

ABFY23.9180225