

Grab Handle/Bed Rail

The Affinity Grab Handle/Bed Rail provides reliable and stable support for users during bed-related activities, including sitting, standing, and transitioning between lying and sitting positions. It anchors securely to the underlying bed frame and removes the potential for shifting during a transfer.



- Support frames slide beneath the mattress and attach securely to the bed base, providing additional stability during use
- Safety straps fasten beneath the bed to prevent shifting and eliminate trip hazards
- Compatible with single and king single bed sizes
- Rail can be installed on either side of the bed
- Easily adjustable in height to accommodate a range of mattress thicknesses, ensuring a secure fit
- Includes a convenient accessory pouch for storing personal items
- Provides sturdy and reliable support, making it easier to get in and out of bed

Grab Handle/Bed Rail BEA013650	Overall Width	Overall Height	Distance between Rails	Strap Length	Product Weight	SWL
	450mm	390 - 520mm	115mm	3100mm	3.7kg	136kg

