SCOOTER SELECTION GUIDE

	Scooter Code	Image		
COMPACT	BOOT SCOOTER (SCT670613)	Exclusive	e	
	BOOT SCOOTER SUPALITE (SCT670616)	Exclusive	e	
	GO-GO® LX + CTS SCOOTER (SCT670612)			
	SONIC SCOOTER (SCT670665)	Exclusive	e	
	QUEST FOLDING COMPACT SCOOTER (SCT670630)			
MEDIUM	GT SCOOTER (SCT670670)	Exclusive	e	
	PATHRIDER® 10 SCOOTER (SCT670615)	\$		
	PATHRIDER® 10 DXE SCOOTER (SCT670507)	Exclusive	e	
	VICTORY 10 LX SCOOTER (SCT670640)	ž: ~~		
	PATHRIDER® 130 XL PURSUIT SCOOTER (SCT670571 / SCT670575)	Exclusive	e	
LAKUE	PATHRIDER® 140 XLE SCOOTER (SCT670505 / SCT670506)	Exclusive	e	
MAXI	RANGERIDER SCOOTER (SCT670581)	Exclusive	e	

Terrain	Battery	Max Range Per Charge	Safe Working Load	Features
Level & smooth terrain. Short distances.	2x 12V. 12AH	11.5 km	118 kg	 Up to 6km/h Disassembles for easy transport
Level & smooth terrain. Short distances.	1x 24V. 12AH	18.5 km	120 kg	 Up to 6.5km/h Disassembles for easy transport Comfort Trac front & rear suspension (CTS)
Level & smooth terrain. Short distances.	2x 12V. 18AH	17 km	136 kg	• Up to 6.4km/h • Disassembles in easy steps
Level & smooth terrain. Short distances.	2 x 12V. 20AH	20 km	136 kg	• Up to 8.8km/h • Great outdoor stability, ease of use and maximum va
For more varied terrain. Short to medium distances.	1x 24V. 11.5AH Li-Polymer	20 km	115 kg	 Up to 6km/h Disassembles for easy transport Easy-to-grip tiller adjustment knob
For more varied terrain. Short to medium distances.	2 x 12V. 35AH	30 km	159 kg	• Up to 10km/h • Front & rear suspension
For more varied terrain. Short to medium distances.	2x 12V. 31AH (40AH optional)	25 km	159 kg	 Up to 8.5km/h Comfort Trac front & rear suspension (CTS)
For more varied terrain. Short to medium distances.	2x 12V. 40AH	39 km	181 kg	 Up to 10km/h Front & rear suspension Pneumatic tyres
Handles rougher terrain well. Long distances.	2x 12V. 40AH	26 km	159 kg	 Up to 9km/h Front & rear suspension Upgrade to 70AH batteries (x2) for improved perform
Handles rougher terrain well. Long distances.	2x 12V. 55AH	25 km	159 kg	• Up to 10km/h • Front & rear suspension • Upgrade to 100AH batteries (x2) for improved perform
Handles rougher terrain well. Long distances.	2x 12V. 100AH	35 km	181 kg	 Up to 10km/h Front & rear suspension High performance and 13" pneumatic tyres
Handles rougher terrain well. Long distances.	2 x 24V. 100AH	30 km	159 kg	 Up to 10km/h Front & rear suspension High performance and 14.5" pneumatic tyres

Features	F		
• Up to 6km/h • Disassembles for easy transport	Aid		
• Up to 6.5km/h • Disassembles for easy transport • Comfort Trac front & rear suspension (CTS)	- lea . p . p		
• Up to 6.4km/h • Disassembles in easy steps	be • Pl re		
• Up to 8.8km/h • Great outdoor stability, ease of use and maximum value			
 Up to 6km/h Disassembles for easy transport Easy-to-grip tiller adjustment knob 			
Up to 10km/hFront & rear suspension			
• Up to 8.5km/h • Comfort Trac front & rear suspension (CTS)	1		
 Up to 10km/h Front & rear suspension Pneumatic tyres 	2 3 4		
• Up to 9km/h • Front & rear suspension • Upgrade to 70AH batteries (x2) for improved performance	5		
Up to 10km/h Front & rear suspension Upgrade to 100AH batteries (x2) for improved performance			
 Up to 10km/h Front & rear suspension High performance and 13" pneumatic tyres 	T U re V		
• Up to 10km/h	S		

our website.



YOUR BEST*

are works closely with Pride Australia to offer a rehensive line of mobility products to help our customers ne most active, independent lifestyles.

- e Scooter products are all registered in the ARTG
- se always read and understand your owners handbook re operating your new mobility scooter
- se take care to observe safety and the speed limit irements of your state



- mfortable front-mounted post seat
- ick-lit battery gauge on monitor
- ont basket as standard
- eek & sporty design
- D headlight, tail light and indicators on some models
- w profile, puncture proof tyres or pneumatic tyres

SAFE USE VIDEO

e are a few safety considerations we'd like you to rstand before operating your device. We strongly mmend all new mobility device customers th our 7 minute YouTube video, **'Getting to**

i a

w your Power Mobility Device'.

can the QR code or visit the 'videos' page on



FLAG

SCA668370

SUN CANOPY

SCA668375

CUP HOLDER

SCA668170

OXYGEN BOTTLE

HOLDER

SCA668330

SCOOTER ACCESSORIES



SADDLEBAG SCA668398



REAR POD SCA668399



REAR BASKET SCA668180



CANE HOLDER SCA668230



DUST COVER SCA668210



MIRROR SCA668300



SCOOTER BAG + L PIECE SCA668340



HOLDER SCA668390

COMPATIBILITY MATRIX

Scan the QR code to view full range of Mobility Scooter Accessories and their compatibility.









MOBILITY SCOOTER **SELECTION GUIDE**







aidacare.com.au

PRIDE[®] MOBILITY SCOOTER RANGE

BOOT SCOOTER

SCT670613



Exclusive

- Feather touch disassembly for easy transportation & storage
- · Front-mounted post seat with arm rests
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller Modular design for easy servicing Colour: Red Length: 1003 mm
- Width: 495 mm Product Weight: 45.6 kg (inc. batteries) Maximum Speed: 6.4 km/h Range: 11.5 km Battery Size: 2 x 12V. 12 AH

SONIC SCOOTER SCT670665



Exclusive

- · Disassembles into 5 light weight pieces
- · Portable and durable scooter
- Puncture-proof tyres
- · LED headlights and brakelights
- Detachable battery with easy access
- Adjustable angle tiller with delta Colour: Red
- Length: 1040 mm
- Width: 495 mm
- Product Weight: 52 kg (incl batteries) Maximum Speed: 6.8 km/h

• Sleek, stylish and easy to operate

- Range: 20 km
- Battery Size: 2 x 12V. 20AH

BOOT SCOOTER SUPALITE SCT670616

- Disassembles into 6 lightweight pieces
- Compact, foldable seat with armrests Easy-to-grip tiller adjustment knob Stable front frame-mounted seat post
- Standard front basket provides storage Colour: White
- Length: 1022 mm
- Width: 480 mm Product Weight: 34 kg
- (incl batteries)
- Maximum Speed: 6.5 km/h
- Range: 18.5 km

Battery Size: 1 x 24V. 12AH Lithium Ion

PRIDE GO-GO[®] LX + CTS SCOOTER SCT670612

1300 133 120 | aidacare.com.au



- Comfort-Trac Suspension
- Lightweight mobility scooter Disassembles into 5 pieces for easy transportation and storage Easy to remove rechargeable battery Soft and easy to adjust tiller Colour: Red Length: 1040 mm Width: 520 mm Product Weight: 39.5 kg
- Maximum Speed: 6.4 km/h Range: 17 km
- Battery Size: 2 x 12V. 18AH

Exclusive PRIDE OUEST FOLDING SCOOTER SCT670630

swL 115kg

GT SCOOTER

SCT670670



- Front and rear LED lighting
- Directional signals
- Deluxe tiller with USB charger
- Cup holder Colour: Red
- Width: 605 mm
- Product Weight: 103 kg (incl batteries)
- Maximum Speed: 10 km/h Range: 30 km
- Battery Size: 2 x 12V. 35AH

PRIDE PATHRIDER® 10 SCOOTER

SCT670615



High performance, sleek and sporty

- Low profile, puncture proof tyres
- Back lit battery gauge
- · Front-mounted post seat with armrests
- Eight-position seat rotation Wrap around easy drive tiller Colour: Blue Length: 1190 mm Width: 565 mm Product Weight: 63.5 kg Maximum Speed: 8.5 km/h Range: 21 km

/ 159kg \

PRIDE PATHRIDER[®] 10 DXE SCT670507

Exclusive

New style front tiller with LED lights

Battery Size: 2 x 12V. 31AH (40AH option)

- Comfortable high back seat USB charging point for devices
- 10" pneumatic tyres
- Integrated storage space with cup holder built inside.
- Colour: Blue
- Length: 1194 mm Width: 560 mm
- Product Weight: 68 kg
- Maximum Speed: 10 km/h
- Range: 39 km
- Battery Size: 2 x 12V. 40AH

separation with only one hand

Colour: Red

Length: 1194 mm

Product Weight: 120 kg

Maximum Speed: 9 km/h

Battery Size: 2 x 12V. 40AH

Width: 565 mm

Range: 26 km

Full LED lighting package includes curb

and rear lights and directional signals

PRIDE VICTORY 10 LX SCOOTER SCT670640



PRIDE RANGERIDER SCT670581

181kg



*Scooter range is manufactured by Pride Mobility and distributed by Aidacare



- Folded & ready for transport in 3 steps · Can be stowed upright or laid down Great performance, indoors and out Articulated front end
 - · Front and rear LED headlights

Product Weight: 23.5 kg

Range: 20 km

Battery Size: 1 x 24V (11.5 AH) Li-Polymer

Exclusive





Advanced suspension for a smooth ride

- Length: 1200 mm

STANDARD POWER WHEELCHAIRS

PATHRIDER® 130 XL PURSUIT

SCT670571 - RED. SCT670575 - YELLOW

- High performance, sleek and sporty
- Large 13" pneumatic tyres
- · LED headlights, brakelights & indicators

· Armrests, head rest & lumbar support Rugged front and rear suspension Colour: Red Length: 1359 mm Width: 635 mm Product Weight: 109 kg (incl batteries) Maximum Speed: 10 km/h Range: 25 km

Battery Size: 2 x 12V. 55AH

WHILL MODEL C2 POWER WHEELCHAIR

PWS655320



 Omni-wheels for tight turning radius Electromagnetic brakes provide safe

- stopping, even on inclines
- Adjustable seat height, back support
- Control remotely with iPhone/Android Anti-sway system
- Colour: Black Lenath: 990 mm Width: 450 mm
- Product Weight: 52 kg
- Maximum Speed: 8 km/h
- Range: 16 km Battery Size: Lithium-ion (25.2V. 10Ah)

PRIDE PATHRIDER[®] 140 XLE SCT670505 - YELLOW, SCT670506 - RED

Exclusive

Exclusive

- New advanced LCD panel display
- High performance ride quality
- Large 13" pneumatic terrain-gripping low-profile tyres
- · Directional signals and large lower headlight for pathway illumination Colour: Yellow or Red Length: 1422 mm Width: 692 mm Product Weight: 125 kg Maximum Speed: 10 km/h Range: 35 km Battery Size: 2x 12V. 100AH

PRIDE LITESTREAM[®] CARBON PWS654185



- Lightweight carbon fibre frame
- 120cm turning radius
- Under seat pouch for extra storage
- Width: 570 mm
- Product Weight: 16 kg (excl battery)
- Maximum Speed: 6 km/h

PRIDE JAZZY SELECT 6 (MID WHEEL DRIVE) PWS654242

- Maximum front stability and in-line motor technology
- 6" front and rear castors
- · Height and depth adjustable foot platform Length & Width: 900 x 600 mm Product Weight: 75 kg Maximum Speed: 6 km/h
- Range: 21 km

PRIDE IGO FOLDING POWER WHEELCHAIR PWS654175



- Quickly folds in 5 easy steps
- Unique seating for added comfort
- Front suspension
- · Large, secure storage compartment
- Easy-to-use VR2 controller Length & Width: 875 x 600 mm Maximum Speed: 5.95 km/h Product Weight: 19.8 kg Seat Width: 470 mm Battery Size: 2 x 6.6 Ah Lithium Ion

PRIDE GO CHAIR - NEXT GENERATION (REAR WHEEL DRIVE) PWS654200RE

- · Disassembling compact design
- Tight turning radius
- 5-point swivel seat with latch release, 52cm lap belt, Length & Width: 838 x 574 mm Maximum Speed: 5.95 km/h Product Weight: Heaviest Part 16.3 kg Seat Width: 460 mm Range: 14 km Battery Size: 2 x 12V. 18AH deep-cycle



To view more mobility aids or complex rehab power chairs, scan the QR code or visit our website



◀ Aspire VOGUE Carbon Fibre Walker



shock absorbers · Reversing sensors, LED light Dual hydraulic brakes Length: 1560 mm Width: 830 mm Maximum Speed: 10 km/h Range: 30 km

• Flip-up height and angle

adjustable armrests



SALES · HIRE · SERVICE





- · Optional manual setting · Folds for easy transport Colour: Black
- Length: 930 mm

- Range: 15 km
- Battery Size: 24V. 12Ah Lithium Ion

- CTAS suspension with adjustable



Seat Width: 460 mm

aidacare.com.au

Battery Size: 2 x 12V. deep-cycle