

KINERIS PRESSURE REDUCING CUSHION RANGE

Ideal for patients at very high risk of pressure injury development. It is suitable for both active and non-active users, accommodating a large range of user heights and weights.



- **Accurate Pressure Adjustment:** Manometer enables exact pressure adjustment every time. Easy for carers!
- **Single-valve:** For simple inflation and maintenance.
- **Reduced Friction and Shearing:** Cell mobility minimises friction and shearing forces.
- **Optimal Height:** Available in low and high profiles.
- **Comfortable and Breathable:** The removable cover is stretchable, and breathable, and aids body temperature regulation during prolonged sitting.
- **Stays in Place:** The inner cover features a non-slip base to keep the cushion secure.
- **Versatile Sizing:** Available in various sizes to accommodate different needs.

WinnCare
Kinéris

Code	Width Depth Height	Profile	Code	Width Depth Height	Profile
PTC627100	16.5in 16.5in 2.75in	Low	PTC627200	16.5in 16.5in 4in	High
PTC627110	16.5in 18.5in 2.75in	Low	PTC627210	16.5in 18.5in 4in	High
PTC627120	18.5in 18.5in 2.75in	Low	PTC627220	18.5in 18.5in 4in	High

Exclusively distributed by