

EQUIPMENT SELECTION MADE EASY




PAEDIATRIC CATALOGUE



# **Paediatric**

- 04 Standing Aids
- 06 Seating Aids
- 10 Strollers
- 12 Walking Aids
- 15 Manual Wheelchairs
- 20 Power-Assist Systems
- 22 Power Wheelchairs
- 26 Head Supports
- Night Time Positioning
- 29 Bathroom



## Standing Aids

#### **Jenx Multistander**

- Versatile standing frame, suitable for children from around 9 months right up to 13 years
- PU boards and supportpads are flexible, easy clean and hygienic
- De-rotational thoracic & pelvis straps/ pads fasten in the centre and adjust independently from either side



- pelvic boards cater for the child's individual proportions
- Deliberately open design ensures you can always see precisely how the user is positioned each and every time



**LECKEY** 

S.W.L

22KG

Code	Description	User Standing Height	Chest Pad Height	S.W.L	Wty
PED479630	Stander - Jenx Standz - Abduction Stander (As Per Script)	1160 - 1750mm	490 - 1170mm	35 - 60kg	2 Years



#### **Squiggles Standing System**

- Squiggles Standing System supports the child in a comfortable, secure, prone, supine or upright standing posture
- Chassis can clip on and off in seconds for folding and storage
- differing needs • Accessories available include tray, knee blocks, supports, headrests which all

• Adjustable to accommodate children with

attach without tools Squiggles Stander, Ages 1 - 5 years 750 - 1100mm 550 - 820mm

NDIS Item Code 05\_053603010\_0103\_1\_2



## **Squiggles+ Stander**

- Hip abduction of up to 60° in prone. upright and supine standing.
- Fully horizontal loading in supine.
- New wrap around knee supports to better accommodate all sizes, with added comfort & support.
- · New, improved, pivot chassis.
- Height-adjustable, foldable push handle for ergonomic use.
- Individually adjustable footplates to accommodate leg length discrepancy.



**LECKEY** 

Code	Description	User Standing Height	Chest Pad Height	Wty
PED479380	Squiggles+ Stander, Ages 1 - 5 years (As per Script)	750 - 1110mm	550 - 820mm	2Years



#### **Mygo Stander**

- The Mygo Stander provides a 3-in-1 prone, supine and upright standing positions
- Comes with height adjustable chest pad that is ergonomically designed to imitate

the shape of the rib cage

 Includes hip support, chest support, 4 castors with brakes and height adjustab footplates and knee supports

	A
le	S.W.L 50-60KG

**LECKEY** 

(0)

S.W.L

50-100KG

**LECKEY** 

Code	Description	Overall Width	Overall Length	Chest Width	Hip Width	Height of User Range	Wty
PED479335*	Mygo Stander, Size 1, Ages 4 - 10 years (As per Script)	630mm	978mm	160 - 280mm	220 - 330mm	1080 - 1450mm	2Yrs
PED479336	Mygo Stander, Size 2, Ages 8 - 14 years (As per Script)	650mm	1130mm	200 - 320mm	280 - 390mm	1350 - 1750mm	2Yrs

1000 - 1400mm

1400 - 1800mm

NDIS Item Code 05\_053603010\_0103\_1\_2

#### **Pronestander**

choice of 5 prone positions ranging from 45° to 85°

**Jenx Monkey Prone Stander** 

introduction to standing, providing prone to

PED479375 Monkey Standing Frame - Jenx NR-50 (As per Script)

· Jenx Monkey is a gentle and fun

• De-rotational thoracic & pelvis straps/

pads fasten in the centre and adjust

Polyurethane support surfaces with

• The Freestander is a robust upright

• Designed to assist in standing and

weight bearing development

stander ideally suited to the classroom,

independently from either side

upright positioning

antimicrobial agent

Freestander

therapy room or home

- belt, angle adjustable tray, activity bowl and a mobility base with 4 locking castors · Assists in standing and weight bearing 470 - 680mm

NDIS Item Code 05\_053603010\_0103\_1\_2



- The Pronestander is a stander with a

Code	Description	User Standing
PED479250*	Pronestander, Ages 1.5 - 5 years (As per Script)	800 - 1150
PED479260	Pronestander, Ages 4 - 10 years (As per Script)	1100 - 1450
PED479270	Pronestander, Ages 9 - 18 years (As per Script)	1400 - 1850

### **LECKEY**

S.W.L

25KG

medifab<sup>®</sup>

· Highly adjustable to suit individual users

• Integral angle-adjustable tray with both

bowl and smooth tray infill for a wide

• Quick, simple & easy to use. Comes as a

complete unit and very few adjustments

variety of different activities

are required

 Include hip belt bowl ar castors

• Highly adjustable to suit individual users

• Includes chest and knee supports, hip



es chest and knee supports,	
t, angle adjustable tray, activity	
nd a mobility base with 4 locking	60
s for stability	

975 - 1350mm

510 - 750mm

Code	Description	User Standing Height	Chest Pad Height	S.W.L	Wty
PED479170*	Freestander, Ages 1.5 - 5 years (As per Script)	800 - 1151mm	470 - 680mm	60 kg	2 Years
PED479180	Freestander, Ages 4 -10 years (As per Script)	1100 - 1448mm	680 - 975mm	80 kg	2 Years
PED479190	Freestander, Ages 9 - 18 years (As per Script)	1400 - 1849mm	975 - 1350mm	90 kg	2 Years

05\_053603010\_0103\_1\_2

**LECKEY** 



2 Years

2 Years



#### **Horizon Stander**

PED479105

PED479130

- Durable 3-in-1 stander which can be positioned in prone, supine or upright
- Fully adjustable to suit the individual user and is ideal in multi user environments

(As per Script)

Horizon, Ages 4 - 10 years

Horizon, Ages 9 - 18 years

Horizon, Ages 14 years - adult

· Available with headrest, adjustable sandals, lateral supports and pelvic positioning belt

820 - 1150mm

1000 - 1450mm

S.W.L	Wty
50 kg	2 Years
80 kg	2 Years
100kg	2 Years



**Standing** 

Aids

# Seating **Aids**





#### **Early Activity System**

**LECKEY** 

**30KG** 

**LECKEY** 

S.W.L

30-50KG

- Strong, durable, practical and adjustable seating system, providing moderate support for kids and young adults aged 1-18 years
- System supplied as a "kit" to support a young child to sit, lie on their tummy, back or side
- Can be used by therapists as well as parents
- Designed to suit children up to 3 Years

Code	Description	Seat Depth	Hip Width	Backrest Height	Wty
PED480100	Early Activity System	190 - 260mm	150 - 200mm	100mm	2 Years



#### **PAL Classroom Seat**

- · Attractive, robust and easy to use classroom and nursery seating system for kids with mild to moderate postural needs
- Matching tables available
- Improves stability, reduces fatigue and longer periods of concentration and fine motor activity

Code	Description	Seat Width	Seat Depth	Backrest Height	S.W.L	Wty
PED478800	PAL Classroom Seat, Size 1, Ages 1 - 3 years (As per Script)	150 - 250mm	180 - 240mm	265mm	30kg	2Yrs
PED478802	PAL Classroom Seat, Size 2, Ages 2 - 5 years (As per Script)	170 - 270mm	220 - 290mm	315mm	35kg	2Yrs
PED478804	PAL Classroom Seat, Size 3, Ages 4 - 8 years (As per Script)	230 - 330mm	250 - 330mm	380mm	40kg	2Yrs
PED478806	PAL Classroom Seat, Size 4, Ages 7 - 12 years (As per Script)	280 - 380mm	300 - 400mm	450mm	50kg	2Yrs



1300 888 052

### **Everyday Activity Seat**

**LECKEY** 

• Designed to help physical, cognitive and sensory development for babies with special needs

- 4-point harness gives a secure, stable base of support
- Height and angle adjustable flip-up footplate to facilitate transfer
- Highly configurable with 3 sizes, a range of chassis options and 5 colour options

Code	Description	Seat Depth	Hip Width	Backrest Height	Handle Height	S.W.L	Wty
PED490010	Everyday Activity Seat - Size 1 (As per Script)	150 - 280mm	160 - 280mm	300 - 400mm	200 - 300mm	30 kg	3Yrs
PED490015	Everyday Activity Seat - Size 2 (As per Script)	230 - 380mm	260 - 380mm	400 - 520mm	200 - 300mm	50 kg	3Yrs

## **Seating Aids**



#### **Squiggles Seating System**

**LECKEY**°

- Squiggles Seating System is designed for younger children requiring postural support
- Adjustable seating system with tilt in space mechanism
- Designed to suit children from 1 to 5 years and allows for growth
- Seating accessories include thoracic supports, sandals, tray, harness plus others
- Can be ordered with a choice of base including wheelchair base, stroller base or high low
- Seating system fits a number of transporting systems







#### **Mygo Seating System**

**LECKEY**°

S.W.L

50KG

- Multifunctional seating system designed for children and adolescents requiring moderate to complex postural support
- Includes leg guides, four-point pelvic positioning harness and adjustable sacral support
- Seat width and depth are continuously adjustable allowing for a large degree of growth
- Available with unique Leckey pelvic cradle for improved postural control
- Range of additional postural supports available

Seating system fits a number of transporting systems								
Code	Description	Seat Width	Seat Depth	Backrest Height	S.W.L	Wty		
PED479339	Mygo S/System, Size 1, Ages 3 - 10 years (As per Script)	200 - 325mm	270 - 420mm	360 - 470mm	50 kg	2Years		
PED479338	Mygo S/System, Size 2, Ages 8 - 14 years (As per Script)	220 - 345mm	350 - 470mm	460 - 570mm	60 kg	2Years		



#### **Phoenix Paediatric Chair**



- Gives maximum support to the upper body through lateral and head supports, also adapt to the changing needs of a person being rehabilitated or meet the long term needs of a progressive condition
- Adjustable lateral, shoulder and head supports can conform to posture & position asily
- Tilt in space and adjustable components. Meets range of individual support needs, sizes, and levels of mobility



**Seating Matters** 


Code	Description	Seat Depth	Hip Width	Backrest Height	Seat to Floor Height	S.W.L	Wty
CHP197800	Phoenix Chair (As per Script)	290 - 380mm	250 - 400mm	680 - 830mm	540mm	160 kg	7 Years (Frame) 1 Year (Upholstery)


NDIS Item Code 05\_180921121\_0103\_1\_2

Available in multiple manual & electric combinations



#### **Sorrento Paediatric Chair**





- Tilt in space, forward tilt and adjustable back angle recline can be used to improve posture and redistribute pressure
- · Optional lateral supports and pillows

Code	Description	Seat Depth	Hip Width	Backrest Height	Seat to Floor Height	Wty
CHP197819	Sorrento Chair	290 - 380mm	250 - 400mm	650mm	540mm	7 Years (Frame) 1 Year (Upholstery)

NDIS Item Code 03\_090000911\_0103\_1\_1





### **Seating Aids**



#### **Snooza Chair**

- All-foam Snooza Chair uses gravity to keep the occupant safe with good postural stability while they are out of their wheelchair
- Soft snug close support helps the occupant to relax & high armrests promote better breathing
- Carry handles make transportation easy & large storage pockets located on each side
- Wipe-clean medical grade vinyl cover.
   The centre comfort wrap supplied is removable and washable
- Standard Snooza Chair kit includes:
- 1 x Wedge for angle adjustment1 x Adjustable headrest
- 1 x Padded hip belt



Code	Description	Floor to Seat	Seat Depth	Back Height	Width At Knees	Width At Hip	Overall Width	Overall Length
PED472000	Snooza Chair, Size 1, Apple	210mm	260mm	610mm	260mm, 320mm (w/infill pads)	200mm, 250mm (w/infill pads)	460mm	610mm
PED472010	Snooza Chair, Size 1, Royal Blue	210mm	260mm	610mm	260mm, 320mm (w/infill pads)	200mm, 250mm (w/infill pads)	460mm	610mm
PED472020	Snooza Chair, Size 1, Orange	210mm	260mm	610mm	260mm, 320mm (w/infill pads)	200mm, 250mm (w/infill pads)	460mm	610mm
PED472030	Snooza Chair, Size 1, Hot Pink	210mm	260mm	610mm	260mm, 320mm (w/infill pads)	200mm, 250mm (w/infill pads)	460mm	610mm
PED472040	Snooza Chair, Size 2, Apple	240mm	320mm	700mm	300mm, 350mm (w/infill pads)	230mm, 280mm (w/infill pads)	480mm	720mm
PED472050	Snooza Chair, Size 2, Royal Blue	240mm	320mm	700mm	300mm, 350mm (w/infill pads)	230mm, 280mm (w/infill pads)	480mm	720mm
PED472060	Snooza Chair, Size 2, Orange	240mm	320mm	700mm	300mm, 350mm (w/infill pads)	230mm, 280mm (w/infill pads)	480mm	720mm
PED472070	Snooza Chair, Size 2, Hot Pink	240mm	320mm	700mm	300mm, 350mm (w/infill pads)	230mm, 280mm (w/infill pads)	480mm	720mm
PED472080	Snooza Chair, Size 3, Apple	270mm	400mm	750mm	370mm, 420mm (w/infill pads)	270mm, 320mm (w/infill pads)	570mm	870mm
PED472090	Snooza Chair, Size 3, Royal Blue	270mm	400mm	750mm	370mm, 420mm (w/infill pads)	270mm, 320mm (w/infill pads)	570mm	870mm
PED472100	Snooza Chair, Size 3, Orange	270mm	400mm	750mm	370mm, 420mm (w/infill pads)	270mm, 320mm (w/infill pads)	570mm	870mm
PED472110	Snooza Chair, Size 3, Hot Pink	270mm	400mm	750mm	370mm, 420mm (w/infill pads)	270mm, 320mm (w/infill pads)	570mm	870mm
PED472120	Snooza Chair, Size 4, Apple	300mm	450mm	810mm	380mm, 430mm (w/infill pads)	300mm, 340mm (w/infill pads)	600mm	910mm
PED472130	Snooza Chair, Size 4, Royal Blue	300mm	450mm	810mm	380mm, 430mm (w/infill pads)	300mm, 340mm (w/infill pads)	600mm	910mm
PED472140	Snooza Chair, Size 4, Orange	300mm	450mm	810mm	380mm, 430mm (w/infill pads)	300mm, 340mm (w/infill pads)	600mm	910mm
PED472150	Snooza Chair, Size 4, Hot Pink	300mm	450mm	810mm	380mm, 430mm (w/infill pads)	300mm, 340mm (w/infill pads)	600mm	910mm



#### Carrot 3000 Car Seat

- Carrot 3000 caters for safety and the importance of keeping the child in a supported posture for travel
- Various extension accessories and m ultiadjustable features allow longer use during the child's life and suit the wide spectrum of different needs
- Ergonomicallydesigned head support holds the head and protects from side impacts and centrifugal forces

### **LECKEY**

- 2 positions of seat depth adjustment for a better fit
- Shoulder protector wings: Extendable wings hold and protect the body from side impacts and centrifugal forces
- Large chest support holds the body securely
- Recline function with angle adjustable seat and backrest

and centini	agai forces
Code	Description
PED478199	Paediatric Car Seat - Carrot 3000 (As per script)



#### **Goto Size 1**

- Adjustable Head Supports Pick a standard headrest for simple contoured comfort or the advanced headrest which offers more sculpted support
- Adjustable Laterals The laterals move high, low, wide or narrow to give a reassuring 'hug', while the head support slides up and down to the ideal height
- Optional Floorsitter -The new and improved Floorsitter has 26 recline positions. It's ideal for finding the exact level of support your little one needs for playtime or a power nap
- Optional Pressure Cushion



Code	Description	Seat base to top of headrest Max / Min	Lateral support width Max / Min	Seat Depth	Age Range
PED490500	GoTo Seat - Advanced Headrest - Light Blue / Dark Blue - Size 1	540mm MAX 380mm MIN	260mm MAX 190mm MIN	175mm	1-3 Years
PED490550	GoTo Seat - Advanced Headrest - Red Vinyl - Size 1	540mm MAX 380mm MIN	260mm MAX 190mm MIN	175mm	1-3 Years

#### permobil

### Goto

aidacare.com.au



#### Goto Size 2

- Adjustable Head Supports Pick a standard headrest for simple contoured comfort or the advanced headrest which offers more sculpted support
- Adjustable Laterals -The laterals move high, low, wide or narrow to give a reassuring 'hug', while the head support slides up and down to the ideal height
- Optional Floorsitter The new and improved Floorsitter has 26 recline positions. It's ideal for finding the exact level of support your little one needs for playtime or a power nap
- Optional Pressure Cushion



Code	Description	Size	Seat base to top of headrest	Lateral support width	Seat Depth	Age Range
PED490600	GoTo Size 2 - Advanced Headrest - Light Blue / Dark Blue with Floorsitter & cushion	Size 2	680mm MAX 510mm MIN	300mm MAX 210mm MIN	244mm	3-8 Years
PED490605	GoTo Size 2 w/Advanced Headre/ Floorsitter/Cushion - Red Vinyl	Size 2	680mm MAX 510mm MIN	300mm MAX 210mm MIN	244mm	3-8 Years
PED490610	GoTo Size 2 - Advanced Headrest - Pink / Purple with Floorsitter & cushion	Size 2	680mm MAX 510mm MIN	300mm MAX 210mm MIN	244mm	3-8 Years

## Scoot 4 in 1



- Scooter Backrest Our improved higher backrest provides kids with more support and reduces the amount of trunk control needed to use Scooot
- Push pole aids the child to propel themselves out of the chair
- Footplates & Wheels Scooot's durable and easily gripped rubber tyres make is suitable for outdoor use. The footplate can be easily adjusted to your child's exact leg length
- Advanced Backrest Cushion The laterals move up and down to give support where it's needed
- Backrest Cushion Lateral Length 790mm MIN 1010mm MAX
- Push Pole Height: H 800mm MIN H 890mm MAX



Code	Description	Backrest width	Backrest height	Seat Base Width	Wheel Arch Width	Hip to end of footplate when seated	Backrest Cushion	Age Range
PED490900	Scoot 4 in 1 with Advanced backrest	260mm	300mm	230mm	230mm	520mm (not extended) 630mm (fully extended)	420mm (Laterals Length)	2-6 Years

### **Playpak Portable Activity Kit**








• Therapeutic and early intervention strategies allow your babies and children to develop healthy habits



Code	Description	Length	Width	Depth
	Folded Bag	440mm	490mm	150mm
	Floor Mat	1340mm	750mm	-
	Long Roll	520mm	90mm	70mm
PED490600	Short Roll	280mm	90mm	70mm
	Positioning Strap	660mm	80mm	-
	Large Horseshoe	440mm	490mm	150mm
	Small Horseshoe	70mm	90mm	200mm







### Strollers

### **Bingo Evo Stroller**

- Modern, Compact folding size and sporty design
- Multifunctional seat unit
- Quickly change the seat from front to rear
- · Seat positioning controls are positioned centrally for easy operation

## medifab

- Individual adjustability thanks to various adjustment options • Perfect individual care thanks to numerous
- Growth capacity in seat width and depth,

back height, and lower leg length



ode	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	
/CA859680	Bingo Evo Cobra Stroller, Size 1, Indoor Frame	180 - 300mm	160 - 340mm	280 - 460mm	160 - 370mm	15kg	1 year



#### **Wonderseat Stroller**

- · Back support shaping allows for accommodation or correction of trunk spinal deformities, that will ultimately be supported through multi-adjustable lateral supports and anterior chest harnesses.
- Adjustable pre-ischial area can be conformed on-the-spot for greater or
- lesser anti-thrust needs
- Strategic positioning pads can be added or removed from underneath the cushion
- Comfortable, padded forearm supports on quick-release hardware
- Feet positioning sandals



spex

	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	
ED491300	Wonderseat - 45° (as per script)	170 - 330mm	180 - 380mm	360 - 480mm	250 - 450mm	2Years



#### **Kimba Neo Stroller**

- Indoor or outdoor stroller for children requiring postural support
- Built-in suspension ensures a smooth ride outdoors
- Easy-to-use, unique folding system for easy transportation or storage
- Includes height-adjustable push handle, angle-adjustable backrest plus tilt-in-space
- A number of seating support accessories available



ottobock.

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Foot- rest	Wty
PED478438	Kimba Neo Stroller, Size 1, Ages 1 - 6 years (As per Script)	200 - 300mm	190 - 310mm	410 - 610mm	190 - 310mm	2 Years
PED478439	Kimba Neo Stroller, Size 2, Ages 4 - 10 years (As per Script)	260 - 350mm	240 - 400mm	570 - 700mm	200 - 370mm	2 Years



aidacare.com.au

#### **WeGo Stroller**

- WeGo has a higher weight capacity than most pushchairs. It carries kids up to 30kg
- It reclines to 140°, so your child can lie back for a little snooze or to be changed and they'll still have the support of the GoTo Seat
- Easily fit your GoTo seat in the seat
- Comes with a full accessory pack, including a sun canopy, rain cover, cup holder, mosquito net and foot muff
- WeGo comes with a full accessory pack, including a sun canopy, rain cover, cup holder, mosquito net and foot muff.
- Latex free
- Outdoor Use
- Breathable Fabrics
- Weighs 11kg

S.W.L 30KG	\

Code	Description	Total Height	Total Weight	Seat Depth	Seat Height	Seat Width	Footrest Shortest Setting	Footrest Longest Setting
PED491200	WeGo - Push Chair with full accessories pack	1120mm	110mm	260mm	500mm	350mm	250mm	330mm

#### **EZ Rider Stroller**

- Ultra lightweight compact folding stroller
- 10 degree seat tilt allows correct positioning for children needing postural support, but also keeps them sitting upright
- · Seating accessories include headrests, torso supports and tray
- Available as crash tested Bus Transport model
- Bus Transport model comes standard with tie-down system, H-Harness, headrest extension, and foot straps
- Available in 4 different seat widths



Convaid

S.W.L

30-77KG

Convaid

Code	Description	Width	Depth	Height	Footrest	Weight	Wty
PED477700	EZ Rider (As per Script)	305 - 455mm	229 - 495mm	451 - 508mm	76 - 559mm	11.4 - 13.2 kg	Limited Lifetime



**Strollers** 

#### **Cruiser Stroller**

- Lightweight compact folding stroller
- Set at 30 degrees seat tilt to allow for more postural support
- · Seating accessories include headrests, torso supports and tray

#### • Available as crash tested Bus Transport model

- Bus Transport model comes standard with
- tie-down system, H-Harness, headrest extension, and foot straps
- Available in 5 different seat widths

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	Wty
PED477260	Cruiser (As per Script)	255 - 455mm	152 - 533mm	533 - 630mm	130 - 580mm	8.6 - 13.6 kg	Limited Lifetime



#### **Rodeo Tilt Stroller**

- Versatile Tilt-In-Space stroller
- Compact folding and lightweight design for easy transporting
- Available as crash tested Bus Transport model
- Available in 4 different seat widths
- Bus Transport model comes standard with tie-down system, H-Harness, headrest extension, and foot straps



Convaid

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	Wty
PED478195	RodeoTilt (As per Script)	255 - 405mm	178 - 457mm	533 - 635mm	230 - 457mm	13.2 - 17.7 kg	Limited Lifetime



## Walking Aids



### **MyWay Walking Aid**

- Open frame design enables children to touch, feel, learn and play at their own pace
- MyWay castors are multifunctional to assist children with their individual needs
- Can be fine tuned to achieve the optimum position for mobility for each child

• Open frame also allows therapists to work closely with the child on their lower limbs for therapeutic benefits

• Frame is foldable, facilitating convenient storage and transportation



**LECKEY**°

Code	Description	User Standing Height	Seat Height	S.W.L	Wty
PED479600	MyWay - Size 1, Ages 1 - 5 years (As per Script)	635 - 920mm	220 - 440mm	25kg	2Yrs
PED479610	MyWay - Size 2, Ages 4 - 11 years (As per Script)	870 - 1190mm	400 - 610mm	50kg	2Yrs
PED479615	MyWay - Size 3, Ages 9 - 16 years (As per Script)	1070 - 1490mm	550 - 810mm	80kg	2Yrs

Walking

Aids

### **Nurmi Neo Walking Aid**

- Innovative walking aid for children
- Positioned behind the child's body whilst walking allowing open directional walking
- Lightweight folding frame is height adjustable to grow with the user

#### ottobock.

THE RID

- Includes reverse roll locks and swivelling front castors
- Options available to customise the frame to the individual user
- Available in 3 sizes



Code	Description	Overall Width	Overall Length	Product Weight	S.W.L	Wty
PED479349	Nurmi Neo Walking Aid - Size 1 (As per Script)	630mm	650mm	5.5kg	25kg	2 Years
PED479350	Nurmi Neo Walking Aid - Size 2 (As per Script)	680mm	760mm	6kg	40kg	2 Years
PED479351	Nurmi Neo Walking Aid - Size 3 (As per Script)	720mm	830mm	6.5kg	55kg	2 Years

Code	Description	Overall Width	Overall Length	Product Weight	S.W.L	Wty
PED479349	Nurmi Neo Walking Aid - Size 1 (As per Script)	630mm	650mm	5.5kg	25kg	2 Years
PED479350	Nurmi Neo Walking Aid - Size 2 (As per Script)	680mm	760mm	6kg	40kg	2 Years
PED//70351	Nurmi Neo Walking Aid - Size 3	720mm	830mm	6 5kg	55kg	2Voare



### **Kidwalk**

- Partially weight bearing, posturally supportive mobility system which encourages mobility and moves with the child
- Easy, tool-free set-up and adjustment
- Child walks unaided while also keeping hands free
- Available in 2 sizes



ottobock.

	Description	Overall Width	Overall Length	Height of User	S.W.L	Wty
ED479355	Kidwalk, Size 1, Ages 1 - 7 Years (As per Script)	640mm	825 - 875mm	610 - 1220mm	30kg	2 Years
ED479356	Kidwalk, Size 2, Ages 5 - 12 Years (As per Script)	700mm	900 - 950mm	1010 - 1520mm	40kg	2 Years

### **Mobilisor (Paediatric)**

• Mobilisor range is a sturdy, folding, height - adjustable walking frame for children

• Includes fixed front wheels and rear glides

Code	Description	Width	Handle Ht.	S.W.L	Wty
WAF704200	Mobilisor Child	440mm	600 - 790mm	60kg	1 Year
WAF708600	Mobilisor Junior	495mm	670 - 960mm	80kg	1 Year



### Nimbo Lightweight Walking Aid

- Lightweight aluminium frame
- Tool free height adjustment
- Soft rubber wheels adhere to any surface and prevent Nimbo from sliding backwards
- One directional rear wheels with aluminium ratchet and pin mechanism
- Latex free handles with flanged ends
- Folds easily for convenient storage and transportation
- Available in 5 sizes Specify size on request
- Accessories available on request



Code	Description	Overall Width	Handle Height	Product Weight	S.W.L	Wty
PED480780	Nimbo Walking Aid -Tyke, Size 1, 390 to 500mm, Goldenrod Yellow	530mm	365 - 470mm	4kg	45.5kg	2 Years
PED480790	Nimbo Walking Aid - Junior, Size 2, 480 to 640mm, Lime Green	555mm	480 - 635mm	5kg	45.5kg	2Years
PED480800	Nimbo Walking Aid - Youth, Size 3, 500 to 770mm, Cornflower Blue	635mm	580 - 770mm	5.5kg	91kg	2 Years
PED480805	Nimbo Walking Aid - Youg Adult, Size 4, 710 to 910mm, Midnight Blue	645mm	720 - 910mm	6kg	91kg	2 Years
PED480807	Nimbo Walking Aid - Adult, Size 5, 920 to 1040mm, Flame Red	740mm	910 - 1040mm	7.72kg	125kg	2Years
PED480810	Nimbo Seat for Lightweight					2Years

### **Walking Tutor Forearm Gutter Frame**

- Tutors the user into an upright gait
- Folding design for storage or transportation
- Height adjustable forearm supports to suit user requirements
- Adjustable hand grips to suit user forearm length
- Four swivel wheels for manoeuvrability
- Powder coated mild steel construction
- Skis or glides available as an option as well as padded chest rest



Code	Description	Depth	Forearm Height	S.W.L	Wty
WAF704400	Walking Tutor - Child	410mm	600 - 790mm	60kg	1 Year
WAF709200	Walking Tutor - Junior	495mm	670 - 960mm	80kg	1 Year
WAF709300	Walking Tutor - Junior w/ Glides	495mm	670 - 960mm	80kg	1 Year
WAF709400	Walking Tutor - Junior w/ Padded Chest Rest	495mm	670 - 960mm	80kg	1 Year

## Walking Aids



#### **Kids Forearm Crutches**

• Just for kids, these light weight yellow Forearm Crutches are height adjustable

Code	Description	Handle Height	Ferrule	S.W.L	Wty
WAC691700	Crutches Forearm Child 4 - 7 years	406 - 559mm	18mm	100kg	1 Year
WAC691600	Crutches Forearm Child 6 - 10 years	457 - 656mm	18mm	100kg	1 Year

03\_120000911\_0105\_1\_1



#### **Underarm Crutches - Red Dot**

- Patented clip lock adjustable handle and push button foot pieces for fast, precise measuring and fitting
- 22mm ferrules

- High density rubber cushions and handgrips providing long wear and greater comfort
- Lightweight aluminium frame, child weighing 1.1kg each, and youth - adult tall weighing 1.25kg each

Code	Description	User Height	Crutch Height	S.W.L	Wty
WAC692300	Underarm Crutches - Child	1220 - 1370mm	787 - 940mm	136kg	1 Year
WAC692500	Underarm Crutches - Youth	1320 - 1570mm	940 - 1140mm	136kg	1 Year



#### **Underarm Crutches - Wing Nut**

- Easy to use Wing-Nut height adjustment
- High density rubber cushions and handgrips
- providing long wear and greater comfort

Code	Description
WAC692503	Underarm Crutches Youth - Small

NDIS Item Code 03\_120000911\_0105\_1\_1 DVA AP03



**aspire** 

### Upsee

sales@aidacare.com.au


- Creates the possibility of social and emotional development by letting children interact with their peers at eye-level
- Upsee can be used to work on a variety of mobility and ability skills, from weight-bearing and walking to improving head control
- When your child is motivated to explore and move and weight-bear, it creates a real possibility for long-term development

· Lightweight aluminium frame

• 22mm ferrules

1370 - 1575mm 940 - 1145mm

Code	Description	Shoulder to Hip	Chest Circumference	Lower Hip Circumference	Max. User Weight	Wty
PED490700	Therapist Kit	-	-	-	-	1 Year
PED490705	Upsee - Small - Blue	310 - 360mm	520 - 600mm	520 - 600mm	15kg	1 Year
PED490710	Upsee - Medium - Blue	350 - 390mm	570 - 650mm	580 - 670mm	20kg	1 Year



# Manual Wheelchairs





### **Aspire Rehab RX JNR Wheelchair**

- 🦪 aspire
- Designed for the increased postural support and pressure care of children, via easy repositioning and advanced adjustability.
- Tilt (26°) and backrest recline (30°) this facilitates weight re-distribution and the preservation of skin integrity.
- Highly adjustable depth and width adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.
- Breathability support surfaces offer multi-layered foam with breathable multi-stretch covers
- · Contoured, built-up backrest and seat cushion the backrest has a gentle lumbar build up and moderate lateral supports; together with the abductor and adductor build-up at the seat cushion, midline positioning is optimised.
- Advanced articulating legrests closely mimic the normal movement at the knee upon elevation

S.W.L 180KG 2 YEAR WARRANTY Frame YEAR WARRANTY

ode	Description	Seat Width	Overall Width	Seat Depth	Seat Height	Backrest Height	Armrest Height	Min. Chair Weight
/WS449630	Rehab RX JNR	320 - 370mm (13"- 15")	600mm	440 - 490mm	Front: 410 - 440mm Rear: 400 - 475mm	460 - 710mm	310 - 460mm	39.5kg



#### **Aspire Evoke 2 JNR Wheelchair**

- The Evoke 2 Junior is designed to promote independent manoeuvring and handling by juniors, while making lifting and handling a breeze for parents and caregivers.
- Swing away and removable arm rests this facilitates user transfers and also makes for easier transport and storage.
- Swing-in / out, height adjustable and removable foot rests - this facilitates user transfers and allows optimal positioning for the changing needs.

- aspire
- cushion and padded backrest. • Depth adjustable seating with rear angle adjustability - for more accurate positioning, customisability and increased user comfort.

• Strong, Nylon upholstery with a comfortable

• Quick release 16" solid rear tyres - the wheels are easily removed to reduce the weight of the chair for transport and storage.



2 YEAR WARRANTY 1 YEAR WARRANTY

			FOL	DING	Com	ponents
Seat Width	Overall Width	Seat Depth	Seat Height	Backrest Height	Armrest Height	Min. Chair Weight
300mm (12")	510mm	300 - 375mm	360 - 410mm	430 - 490mm	205 - 225mm	14kg

MWS449560 Evoke 2 JNR

NDIS Item Code 05\_122203107\_0105\_1\_2



### **Manual Wheelchairs**



#### Zippie X'Cape Wheelchair

- Rigid wheelchair performance and folding wheelchair convenience using the New Xlock
- Infinite footrest adjustment using the new Z-Finity system



· Versarail for Limitless seating and accessory options



Code	Description	Seat Width	Seat Depth	SeatTo Floor Height	Unit Weight
MWC427501	Zippie X'Cape	203 - 406mm	203 - 457mm	330 - 482mm	8.6kg



#### **Quickie 2 & Quickie 2 Lite Wheelchairs**

OUICKIE

- Folding ultra lightweight manual wheelchair featuring a modular design
- · Versatile, Modular and lightweight
- Over 64,000 possible rehab configurations
- TruFit program



Code	Description	Weight	Seat Width	Seat Depth	Floor-to-Seat Height
MWC427230	Quickie 2 (As per script)	13kg	305 - 508mm	305 - 457mm	360 - 560mm
MWC427235	Q2 Lite (As per script)	8kg	356 - 508mm	356 - 508mm	360 - 560mm

**DVA** AP19



### **Avantgarde Teen 2 Wheelchair**

Aluminum side panels - Light and sturdy

• Knee lever brake

- adolescents and small adults featuring a folding, aluminium alloy frame and aluminium side panels correct width and height for the attendant
- Excellent handling characteristics and is easy to drive with little force
- · Abduction frame creates a good starting position for a correct pelvic position
- A lightweight, self propelled wheelchair for children • Curved push handles can easily be adjusted to the



ottobock.

	707
	S.W.I
NG	90K0

FOLDI

Code	Description	Seat Width	Seat Depth	Weight
MWC427115	Avantgarde Teen (As per script)	220 - 360mm	240 - 400mm	10 kg



16 1300 888 052

#### **Start Junior Wheelchair**

- The first lightweight children's wheelchair that features the same models and accessories as the adult variants.
- Full range of accessories to match the needs of the child
- Lightweight paediatric chair

#### • Uniform caster forks for all wheel sizes

- Durable and resistant upholstery
- Narrow folded dimension with folding lock
- Footrest for small lower leg length
- Available in a range of colours



ottobock.

Code	Description	Seat Width	Seat Depth	Floor-to-Seat Height	Weight
MWC427410	Start Junior 30cm (As per script)	280 - 380mm	240 - 400mm	395 - 520mm	13.5kg
MWC427420	Start Junior 35cm (As per script)	280 - 380mm	240 - 400mm	395 - 520mm	13.5kg

### **TiLite Pilot**

## permobil

Ultra-lightweight design

Carbon fiber frame panels make it 50% lighter and a breeze for kids to navigate

Built-in forward growth adjustment

To grow with your child and maintain optimal rear wheel access

Clean frame design

Features a seat pan and side panels that provide easy access to the adjustment bolts, making changes simple and intuitive



Acid Green Ocean Blue Metallic Metallic

Sapphire Black

Red Pepper





Code	Description	Seat Width	Seat Depth	Floor-to-Seat Height	Weight
MWC427483	TiLite Pilot (As per script)	203 - 355mm (8 - 14")	228 - 381mm (9 - 15")	76 - 368mm (3 – 14.5")	4.9kg



#### **Little Wave XP Wheelchair**

- Ultra-light rigid aluminium wheelchair
- Crescent CG system keeps the side frame free of excess hardware while providing infinite adjustability
- Ki\_Mobility
- Tubular component system provides individual seat height settings



Code	Description	Seat Width	Seat Depth	Seat Height
MWC427039	Little Wave XP (As per Script)	250 - 480mm (10-19")	250 - 410mm (10-16")	290 - 510mm (11.5 - 20")



#### **Little Wave CLIK Wheelchair**

Ki\_ Mobility



- Built in adjustment
- 5th wheel option for extra balance
- Swing away arms

75KG	1
ackrest Height	

Code	Description	Seat Width	Seat Depth	Seat Height	Backrest Height
MWC427493	Little Wave CLIK (As per Script)	200 - 460mm (8-18")	200 - 460mm (8-18")	290 - 510mm (11.5-20")	229 - 457mm (9-18")

### Manual Wheelchairs



#### **Rogue XP**

- Provides maximum adjustment and growth ability so fit, function and ride are optimized.
- Integrated Swing-Away Armrests move with the backrest when adjusting or growing
- 1 3/8" 7000 series aluminum to provide a durable ultra-lightweight solution that meets the active lifestyle needs



adjust both the range and spring rate

• Spinergy Backrest Release Cable

• Maxx Performance is our new spoke wheel that offers performance and style without a high price tag



707
S.W.L 75KG

Code	Description	Seat Width	Seat Depth	Front Seat Height	Rear Seat Height
MWC427048	Rogue XP (As per script)	250 - 460mm	250 - 510mm	340 - 520mm	290 - 510mm



#### **Zippie Zone Wheelchair**

- · Ultra Lightweight, option rich Paediatric rigid chair • Built in Growth Kit and Truefit growth Program



zippie.

Code	Description	Seat Width	Seat Depth	Seat To Floor Height	Unit Weight
MWC427487	Zippie Zone (As per Script)	200 - 460mm	200 - 480mm	360 - 480mm	6.5kg



#### **TiLite Twist**

- Tru-Fit System opens up the world of adjustability by incorporating an infinite adjustment mechanism
- Adjustable Camber Mount incorporates TiLite's Tru-Fit system with up to 3.5" of rear seat height

<ul> <li>Speedloader Adjustable Front Caster Housing</li> </ul>	J
quickly and accurately adjusts fork angle	

Available in a wide range of colours



zippie.

permobil

Code	Description	Seat Width	Seat Depth	Floor-to-Seat Height	Unit Weight
MWC427472	TiLiteTwist (As per script)	203 - 381mm	203 - 406mm	76 - 368mm	4.9kg

### **Zippie Iris Wheelchair**

Rotation-in-space wheelchair

Offers a superior combination of positioning, maneuverability, and increased portability with the XLOCK® folding crossbrace option

sales@aidacare.com.au

This tilting wheelchair uses patented rotation-in-space technology to rotate the seat frame around the user's center of gravity

18 1300 888 052



Shortest Possible Wheelbase

This rotational system provides the shortest possible wheelbase without compromising stability to ensure a smooth and controlled rotation


> Comes with the True Fit Program A complete growth package for free. Cutting-edge rotation-in-space technology designed just for kids



Code	Description	Seat Width	Seat Depth	Floor-to-Seat Height	Weight
MWC427500	Zippie Iris (As per Script)	254 - 457mm	254 - 457mm	320 - 500mm	17 kg



#### **Little Wave Arc**





 Attendant Wheel Lock - Easy access to the wheel lock for the attendant



- Tilt options Tilt the occupant up to 45 degrees.
- Folding recline back
- Folding Adjustable Handle
- Folding Height Adjustable Rigidizer Bar

7 7
S.W.L 75KG

Ki\_Mobility

Code	Description	Seat Width	Seat Depth	Seat To Floor Height	Unit Weight
MWC427015	Arc Little Wave Wheelchair (As per script)	250 - 460mm	300 - 510mm	370 - 480mm	10kg



#### **Little Wave Flip XP Wheelchair**

- Pediatric tilt-in-space solution designed to look more like a stroller with the technology of a wheelchair
- Over 2,000 color combinations

Ki\_ Mobility • Simple "flip" and fold like a stroller, Little Wave

Flip transforms into a small package that can fi in the trunk of the tiniest car.



ode	Description	Seat Width	Seat Depth	Seat Height	Tilt Range
IWC427553	Little Wave Flip XP (As per script)	250 - 460mm (10-18")	300 - 510mm (12-20")	340 - 510mm (13.5-20")	0-45°



#### **Freedom NXT**

- The NXT is the most rigid, folding, paediatric tilt- When the growth spurt hits and the frame in-space manual wheelchair on the market. It is strong and durable, yet compact and lightweight
- The NXT ensure both the child and carer remain comfortable by providing a smoother ride and a user friendly solution

won't fit, the NXT adjusts easily for a good fit today and tomorrow



MWC449180 Freedom NXT (As per script) 255 - 510mm 305 - 510mm





# **Power-Assist Systems**





#### **Smoov One by Alber**

- · Unobtrusive, high quality design
- Robust, rustproof housing and rustproof fork
- Integrated lithium ion battery pack for a range up to 20km, and speed of uyp to 10km/h

## **VEOMS**

- Ergonimic carry handle
- USB-C socket to charge the control unit or

Code	Description	Overall Length	Overall Width	Overall Height	Product Weight
MWC449500	Smoov One by Alber	385mm	145mm	340mm	7.2kg



#### Alber E-Fix E35 & E36

- Lightweight and compact add-on drive can be quickly removed from the wheelchair without the need for tools
- The modern, flat design of the control unit offers simple operation and perfect ergonomics.

- Range of up to 16km with 6.0 Ah battery pack
- Weight of the person of up to 120kg
- Electronic automatic charge and switch off mechanism

Code	Description	SWL
MWC449520	Alber E-Fix E35	120kg
MWC449525	Alber E-Fix E36	160kg



#### **Alber Twion T24**

sales@aidacare.com.au

- Easily transported in a car or on any public service
- Motor, electronics and battery pack integrated in the wheel hub, 24" Wheels
- Stub axle button made from stainless steel
- Patented EasyConnex charging plug uses magnetic force to automatically find the right plug position
- Lithium-ion batteries provide greater ranges and more charge cycles (up to 1,200)
- Small and handy battery charger
- · Charges both batteries similtaneously, charge and switch-off mechanism
- Compatible with Twion mobility app



Negistered Provider

alber

Code	Description	Speed	Range	Product Weight
MWC449510	Alber Twion T24	10km/h	20km	12.6kg

## Power-Assist Systems

#### **Light Drive<sup>2</sup>**

- Instant installation
- Light and easily transportable
- Adapts to all wheelchairs
- Can be transported by plane (12Ah Lithium Battery)
- Removable

- Discreet and unobtrusive design
- · Colour of your choice
- Tackles almost any terrain
- Indoors/outdoors
- Max. user weight: 120kg



Code	Description	Charging Time	Range	Product Weight
WCA858507	Light Drive <sup>2</sup> - Folding	5 - 6 hours	14 - 20km	4kg
WCA858512	Light Drive <sup>2</sup> - Rigid	5 - 6 hours	14 - 20km	4kg



#### **Power Assist**

- Easily attaches directly on to manual wheelchairs
- Features an easy-to-use squeeze trigger mechanism
- Customisable speed settings

- **PRIDE**
- Designed to work in conjuction with most wheelchairs with seat widths of 406 - 508mm
- Maximum speed of up to 5km/h



Code	Description	Suits
VCA858501	Power Assist	Seat Widths 16" - 20"

#### **SmartDrive MX2+ Push Tracker**


### permobil



- Provides great freedom and power to manual chair users
- Ability to travel up hills and over thick carpet
- Lightweight and easy to hook on to manual wheelchairs
- PushTracker app available for iOS and Android devices
- Folding wheelchair compatibility simply pull back on the spring-loaded flanges to quickly attach or remove your SmartDrive MX2
- PushTracker monitor with wristband allows easy tracking of day to day propulsion metrics - including push count, coast time, and more
- Includes battery and mounting hardware

Code	Description	Overall Length	Overall Width	Overall Height	Product Weight
WCA858508	Smart Drive MX2+ Push Tracker	389mm	141mm	242mm	5.7kg





- Selectable speeds for easy transitions between
- Easy to use ergonomic joystick controller
- Full function attendant control option with innovative
- glide to steer function



- Brush-less direct drive wheels add power to our manual wheelchairs
- The Powered Wheels are removable and can be replaced with standard wheels

Code	Description	Charging Time	Range	Product Weight
WCA858506	iXpress Power Add On (As per script)	3 hours	16km/h	17.4kg

DVA AP25





# Power Wheelchairs



### **Explorer Mini**

- The Explorer Mini is a regulatory cleared powered mobility solution that facilitates self-initiated movement and early exploration for young children 12 – 36 months of age with mobility impairments.
- Intuitive design featuring integrated table and backrest is supportive and allows freedom of movement. While the adjustable saddle seat supports hip development and facilities functional sitting or can be removed for incremental standing activities.

## perm<sub>o</sub>bil

- Designed to support development, the Explorer Mini is an ergonomic device with multiple weight bearing surfaces to help promote safe, stable upright postures while providing opportunities to improve strength, endurance and postural control.
- Midline proportional joystick with bright yellow ball handle helps promote visual development, reaching, grasp and release

Code	Description	Max. User Height	Total Chair Height	Seat to Floor Height	Ground Clearance	Base W x L	Max Speed
PED491500	Explorer Mini (As per Script)	1000mm	730 - 930mm	307 - 550mm	250mm	490 x 635mm	2.4 kph

#### **K450 Power Wheelchair**

- Unique 'terra' function allows user access to the ground
- Advanced suspension system



permobil

<ul> <li>Fully growable seating sy</li> </ul>	stem
---	------

Code	Description	Seat Width	Seat Depth	Seat To Floor Height	Base W x L	Max Speed
PWC652554	K450 (As per Script)	250 - 410mm	250 - 460mm	360 - 460mm	610 x 890mm	8 kph

#### **Koala Power Wheelchair**

- True paediatric size power base and seating Tilt
- Fully growable seating systems

•	Tilt	and	elevate	standard



permobil

Code	Description	Seat Width	Seat Depth	SeatTo Floor Height		Max Speed
PWC652563	Koala (As per Script)	250 - 360mm	250 - 360mm	375mm	585 x 810mm	6.5 kpł

### **Power Wheelchairs**

#### **Q300 Power Chair**

- Ultra-compact mid or rear wheel drive bases for indoor agility and outdoor stability
- Compatible with Sedeo Pro & Sedeo Lite seats
- Entry level chairs



SUNRISE MEDICAL.

Code	Description	Seat Width	Seat Depth	Max Speed
PWC655255	Q300 Power Chair, MWD Base (As per Script)	385 - 510mm	410 - 560mm	10kph
PWC655260	Q300 Power Chair, R Wheel, Sedeo Sling Back, VR2 Controller (As per Script)	385 - 510mm	410 - 560mm	10kph
PWC655270	Q300 Power Chair, R Wheel Tilt, Sedeo Padded Back, VR2 Controller (As per Script)	380 - 560mm	400 - 560mm	10kph
PWC655280	Q300 Power Chair, M Wheel, Sedeo Sling Back, VR2 Controller (As per Script)	385 - 510mm	410 - 560mm	10kph
PWC655290	Q300 Power Chair, M Wheel Tilt, Sedeo Padded Back, VR2 Controller (As per Script)	380 - 560mm	400 - 560mm	10kph
PWC655295	Q300 Power Chair, Captain Seat, VR2 Controller (As per Script)	460mm	460 - 510mm	10kph



#### **Q500 Mini Power Chair**

- Available in Mid, Forward, Rear & Hybrid Wheel Drive
- Choice of rehab seating systems provide support and positioning to users from mild to medium clinical needs
- Robust outdoor chair with advanced electronics
- Compatible with Sedeo Pro & Sedeo Lite seats



SUNRISE MEDICAL.

Code	Description	Total Width	Total Length	Seat to Floor Height	Max Speed
PWC655330	Q500, Rear Wheel Drive (As per Script)	640mm	1180mm	415 - 490mm	10kph
PWC655335	Q500, MID Wheel Drive Sedeo Lite	610 - 620mm	1280mm	415 - 490mm	10kph
PWC655340	Q500, MID Wheel Drive Sedeo Lite (As per Script)	610 - 620mm	1280mm	415 - 490mm	10kph
PWC655350	Q500, Front Wheel Drive (As per Script)	640mm	1340mm	415 - 490mm	10kph



### Skippi

#### • The Skippi is a fun, exciting, paediatric powerchair.

- While our standard seat can be made as small as 24cm, it can also be paired with the much loved Squiggles and Mygo seating systems from LECKEY.
- The Skippi Plus model allows for specialty controls, ensuring even the most complex clients' needs can be met



ottobock.

ode	Description	Seat Width	Seat Depth	SeatTo Floor Height	Base W x L	Max Speed
WC652265	Skippi (As per Script)	300 - 380mm	280 - 380mm	430mm	570 x 850mm	6kph



### **Power Wheelchairs**



#### **K300 PS JNR**

- Compact base
- Adjustable suspension
- R-net electronics
- Fully growable seating system



- Height and angle adjustable
- Flip back armrests • Removable leg rests



Code	Description	Seat Width	Seat Depth	SeatTo Floor Height	Base W x L	Max Speed
PWC652500	K300 PS JNR (As per Script)	280 - 410mm	250 - 460mm	420mm	610 x 1020mm	10kph



#### **M300 PS JNR**

- Compact base
- Climbing and traction (CAT) link suspension
- R-net electronics
- Fully growable seating system

perm <sub>o</sub> bil

- Height and angle adjustable
- Flip back armrests
- Removable leg rests



Code	Description	Seat Width	Seat Depth	SeatTo Floor Height	Base W x L	Max Speed
PWC652550	M300 PS JNR (As per Script)	280 - 410mm	250 - 460mm	430mm	610 x 925mm	10kph



#### **Q6 Edge 3.0**

- Upgraded SRS (Smooth Ride Suspension)
- 4.5 mph with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Made in the USA of foreign and domestic parts



QUANTUM

Code	Description	Seat To Floor Height	Max Speed
PWC652223	Q500 Power Chair (As per Script)	419 - 495mm	7.24kph

## **Edge 3 Stretto**

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- Toolless removal of rear shroud with battery access door
- TRU-Balance® 3 seating configurations
- Available with Accu-Trac Advanced Tracking Technology



QUANTUM

Code	Description	Seat To Floor Height	Max Speed
PWC652238	Edge 3 Stretto (As Per Script)	406 - 464mm	6.25kph



#### 4 Front

- CASE (Castor Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Q-Logic 3 EX Advanced Drive Controls
- Advanced stability for smooth obstacle transitions
- Quiet ride technology

ode	Description	SeatTo Floor Height	Max Speed
NC652222	4 Front Power Chair (As per Script)	387 - 500mm	9.6kph



#### **V6 Frontier**

- Greatest off-road capability in the Frontier V6 range
- Exceptionally smooth ride due to the low pressure tyres
- Best Frontier V6 model for kerb climbing
- Splash resistant (not waterproof)



Code	Description	Seat To Floor Height	Max Speed
PWC655050	V6 Frontier (As per Script)	435mm	10kph



#### TDX SP 2

- TDX SP (Total Driving Xperience) takes all the best features of the TDX series and improves upon them
- Designed to meet the needs of rehab clients who demand performance and style
- Enhanced performance
- · Superior ride quality

- Built-in technologies such as a True Centre Wheel drive, Enhanced SureStep, Traction Control Design, Quiet Stability Lock, MK6i Electronics
- Powerful 4-pole motors
- Quieter chassis




Code	Description	Seat Width	Seat Depth	Seat To Floor Height	Base W x L	Max Speed
PWC652317	TDX SP 2 (As per Script)	305 - 558mm	305 - 609mm	438, 463, 489mm	901 x 609/647mm	10/12kph



aidacare.com.au

#### **3G Torque**



• It offers many of the features and options found on the other Storm Series models while maintaining a truly affordable price



Code	Description	Seat Width	Seat Depth	SeatTo Floor Height	Base W x L	Max Speed
PWC652415	3GTorque (As per Script)	300 - 550mm	300 - 550mm	420 - 530mm	630 x 800mm	9kph





QUANTUM



**Night Time Positioning** 





### **Hensinger Head Support**

- Hensinger features a soft foam wrap around the collar to support the head by using base support areas at the occiput and jaw line
- Hensinger head support must be used in conjunction with a chest support
- Hensinger enhances more social interaction and awareness during daily activities
- This head support promotes proper head

Code	Description	Inner Circumference
FBS279900	Hensinger head support - X-Small	250mm
FBS279901	Hensinger head support - Small	300mm
FBS279902	Hensinger head support - Medium	350mm
FBS279903	Hensinger head support - Large	400mm

NDIS Item Code 05\_060315111\_0135\_1\_2

positioning resulting in improved breathing, swallowing and simplified feeding



#### **Sleepform**

- Sleepform has been designed to ensure the comfort of users at night or during rest by inhibiting abnormal postural reflexes, and preventing heat and pressure build ups
- Strong, durable, practical and adjustable sleeping system for kids aged 0-18 years

- **LECKEY**°
- Suits all sleeping positions.
- Grid sheet to record best sleeping positioning • Mouldable sleeping system for cots and
- Temperature control sheet available

single beds

Code	Description
PED477010	Sleep Positioning System - Kit Size 1
PED477020	Sleep Positioning System - Kit Size 2
PED477030	Sleep Positioning System - Kit Size 3
PED477040	Sleep Positioning System - Kit Size 4

### **Paediatric Position Pillows**







For more Paediatric position pillows, visit www.aidacare.com.au to check out the complete Snooooooze range.



### **Simple Stuff Works Packages**



• The Simple Stuff Works way provides high quality education and therapeutic positioning equipment with an aim to eradicate the preventable problems associated with body shape distortion

Code	Description
PTA534400	Cushion - Positioning - Simple Stuff Work (As per script)

#### **Care Wave Cushions**

### medifab<sup>®</sup>



Code	Description	One Size
PTA535075	Cylindrical	700 x 170mm
PTA535095	Semi-Fowler	700 x 500 x 140mm
PTA535080	Decubitus	300 x 1800mm
PTA535090	Heel	900 x 180mm
PTA535095	Delta	550 x 300 x 300mm









#### **Norbed Kid**

- Foldable doors have a top latch lock, which is easy for the carer to operate
- Designed to ensure there is no risk of injury to the child while in bed • Handset buttons are colour coded for ease of
- function recognition • The magnet safety locking mechanism
- ensures correct and safe use by only the carer • The handset can be easily, and safely stowed in its holder and locked at the base of the bed,
- to ensure only the carer operates the bed • From as low as 30cm for self-transfer for able bodied users, to 80cm carer working height






- Mattress platform length designed to anthropometric measurement, for increased support and exceptional comfort of the child
- The mattress platform can be configured with no profiling functions, combined manual rastofix and electric handset functionality, or fully electric operated - to provide optimum support for the child





Description	Width Outer/Inner	Length Outer/Inner	Height adjustment	Under bed clearance	Backrest angle	Thigh angle	Legrest lift	Product weight (without bed ends)	SWL
Norbed Kid - Small	900/800mm	1400/1460mm	300 - 800mm	150mm	0 - 70°	0 - 12°	0 - 24°	86kg	70kg
Norbed Kid - Medium	900/800mm	1600/1660mm	300 - 800mm	150mm	0 - 70°	0 - 12°	0 - 30°	95kg	70kg
Norbed Kid - Large	900/800mm	1790/1850mm	300 - 800mm	150mm	0- 70°	0 - 12°	0 - 30°	97kg	70kg

## **Bathroom**





### **Firefly GottaGo**

- GottaGo is an innovative, compact, portable toilet seat which provides you with the consistency needed for toileting success
- Complete with a neat backpack for stress-free travel and a toilet training guide, making toilet time easier for everyone
- Retractable legs for easy folding, stowing and different height settings

•	Can	be	used	over	any	toilet	or	potty	

- Snug hug-like fit and warm touch materials for extra comfort
- Accessories include the footplate, straps, seat strap, potty and backpack
- Lightweight and sturdy

Code	Description	Seat Width	Seat Depth*	Back Height
PED491155	GottaGo Size 1 w/Seat - Frame and all Accessories	225mm	220mm	240mm
PED491160	GottaGo Size 2 w/Seat - Frame and all Accessories	270mm	290mm	360mm
PED491165	Toilet Seat Accessory - GottaGo Size 1			
PED491170	Toilet Seat Accessory - GottaGo Size 3			

<sup>\*</sup> Minimum Seat Depth



aidacare.com.au

#### **Splashy**

- Felxible use in an environment
- Choose from 26 recline options, from upright at 106 degrees to a fully reclined at 140 degrees. The floorsitter can also be used with your GoTo Seat
- The Splashy harness adapts to your child's postural support needs, with the option of a 5-point or the 3-point pelvic support
- Splashy's warm-touch, quick-dry foam covering feels soft whether it's wet or dry
- The sculpted shape provides a comfy and supportive sitting position
- The harness and bumper supports are easily attached using toggles, meaning you can adjust Splashy's supports to fit your child as they continue to grow
- 4 x detachable 'bumper' supports for use anywhere on the seat where support is required (i.e. head rest, pommel or anti-slide shelf)
- Age range: 1-8 years

S.W.L 30KG	
loorsitter	

**(0)** 

S.W.L

30KG

Code	Description	Seat Back Rest		Seat Base		Floorsitter	
oue -		Height	Width	Depth		Width	Height
PED491000	Splashy - Blue Base with Green Harness	610mm	385mm	300mm	350mm	420mm	630mm
PED491050	Splashy - Blue Base w/Blue Harness	610mm	385mm	300mm	350mm	420mm	630mm



### Bathroom



#### **Splashy Big**

- Felxible use in an environment
- Choose from 26 recline options, from upright at 106 degrees to a fully reclined at 140 degrees
- The Splashy Big harness adapts to your child's postural support needs, with the option of a 5-point or the 3-point pelvic support
- 4 x detachable 'bumper' supports for use anywhere on the seat where support is required (i.e. head rest, pommel or anti-slide shelf)



• Age range: 8-16 years



Code	Description	Seat Back Rest		Seat Base		Floorsitter	
		Height	Width	Depth	Width	Width	Height
PED491100	Splashy Big - Blue Base with Coral Harness	720mm	450mm	400mm	450mm	420mm	795mm
PED491150	Splashy Big - Blue Base w/ Lime Harness	720mm	450mm	400mm	450mm	420mm	795mm



#### **Advance Bath Chair**

- Supports and holds the child to assist with bathing
- Easy-to-use adjustable legs can be set up prior to use
- Simple lever can be used to adjust the backrest recline position in the bath

	LEC	KEY

- Folds flat for easy storage and transportation
- Washable, mesh fabric is breathable and soft to touch
- Accessories include head support, leg supports and shower trolley

A	
S.W.L 72KG	

Code	Description	Seat Width	Seat Depth	Back Height	Warranty
PED478480	Advance Bath Chair Size 1, Ages 1 - 5 years (As per Script)	425mm	225mm	350mm	2 Years
PED478500	Advance Bath Chair Size 2, Ages 4 - 9 years (As per Script)	425mm	275mm	400mm	2 Years
PED478510	Advance Bath Chair Size 3, Ages 8 - 14 years (As per Script)	425mm	335mm	460mm	2 Years
PED478530	Advance Bath Chair Size 4, Ages 12 - 18 years (As per Script)	425mm	400mm	550mm	2 Years



#### **Otter Bathing System**

- Durable plastic frame is lightweight for carrying with easy to clean surfaces
- Soft fabric gives comfortable support and is easily removable for machine washing
- Seat fabric can be slipped off to bathe a child sitting on bottom of bath, supported by the chair back
- Adjustable slip resistant legs
- Height adjustable
- Leg straps enable abduction and adduction for optimum postural control

- Two positioning straps secure child for safety
- Head cushions are height and width adjustable without tools and provide support and comfort
- Can also be used as lateral supports for maintaining straight-line upper body balance
- Folds flat for easy storage and transportation
- Available in 3 sizes Specify size on request
- · Accessories available on request



Code	Description	Seat Width	Seat Depth	Back Height	Warranty
PED480829	Otter Bathing System, Small	370mm	330mm	450mm	2 Years
PED480830	Otter Bathing System, Medium	370mm	330mm	640mm	2 Years



#### **Sharky Bathroom Chair**

- Versatile, practical, and functional! These aspects have been the primary focus throughout the development
- Backrest recline
- 45° Tilt in Space
- Hygiene Opening Removeable section in rear area to assist with hygiene
- Height Adjustable
- Compact Folding Folds down to a height of 28cm for easy transport and storage
- Removable Padding for easy care and cleaning



Code	Description	Overall Width	Overall Depth	Back Height	Seat Angle	Backrest Angle	Knee Angle	Footplate Angle	Wheel	Wty
BTS105970	Sharky Bathroom Chair (As per script)	240 - 350mm	270 - 400mm	420mm	-10° to 45°	72° - 120°	78° - 138°	80° - 100°	75 - 100mm	2 Years

#### **RAZ - CAT Mobile Shower Commode**



- For pediatrics, adolescents or smaller adults, the compact Raz-CAT provides a fit that is better suited to physically smaller clients
- The base and seat frames retain the 18" width of the Raz-AT so that the chair can still be positioned over a toilet
- The seat frame depth has been reduced and a depth-adjustable back offers depths from 12" to 15" (to front of frame) or up to 17" with the front of the seat overhanging the seat frame
- Custom seats are available

- The standard seat has an elongated aperture that works with seat depths from 12" to 15"
- The base frame has been shorted by 3", which makes the chair more maneuverable
- The depth-adjustable back provides growth capability essential component for pediatrics
- Back frames can be ordered in widths 14", 16" and 18"
- Available back heights are 15" and 18"



Code	Description	Overall Width	Overall Length	Seat Depth	Back Height
BTC091800	RAZ - CAT Mobile Shower Commode (As per Script)	550 - 460mm	810mm	280 - 400mm	380 or 460mm









