

# Permobil

## M1



**Contact Aidacare for  
more information on  
1300 133 120**



### **The M1 Core**

**Exclusive to Aidacare from RRP \$9750\***

#### **One chair. Endless Options.**

Finally, a true Permobil quality solution for those users who require a more traditional rehab style seating system with limited to no power seat functions. Designed with adaptable seating and positioning needs in mind, the M1 offers all the flexibility without the compromise.

#### **Proven Technology.**

Built on the modular chassis of the Permobil M3 Corpus®, the M1 features the same reliable drive units, patented FlexLink™ suspension and precision engineered CasterGlide™ technology to help ensure smooth, confident operation in a variety of driving scenarios.

#### **VR2 System.**

The VR2 combines ergonomic elegance with a stylish design, offering straightforward and intuitive controls for users. Designed with user convenience in mind, this system ensures easy maintenance and serviceability, keeping users at the forefront of comfort and control.

\*Indicative model and pricing

Excludes dealer delivery, associated costs, and any additional accessories required for your needs.  
Contact Aidacare for complete pricing and options.

**permobil**

# Tilt at the push of a button.

The power tilt option allows tilting the seat up to 50° to promote a resting position or pressure distribution. Tilt can also be used to improve comfort, trunk stability and head control.



## Specifications

<b>Maximum user weight</b>	136 kg <sup>††</sup>
<b>Maximum speed</b>	8 km/h
<b>Range</b>	25 km*
<b>Base width</b>	24"
<b>Base length</b>	35.5"
<b>Minimum turning radius</b>	20.5"
<b>Ground clearance</b>	3"
<b>Recommended battery type</b>	Group 34
<b>Weight including batteries</b>	145 kg <sup>**</sup>
<b>Drive electronics</b>	VR2 70A (Advanced)
<b>Seat-to-floor height</b>	16"

<b>Tilt options</b>	0° – 50° posterior
<b>Recline options</b>	85° – 110° (manual)
<b>Legrest angle</b>	80° – 135° (manual)
<b>Backrest height</b>	18" – 29" (by 1" incr.)
<b>Seat depth</b>	14" – 20" (by 1" incr.)
<b>Seat width</b>	14" – 20" (by 2" incr.)
<b>Armpad options</b>	4"x13"

<sup>†</sup>Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.  
<sup>††</sup>The maximum user weight for a seat depth of 14" cannot exceed 75 kg.

\* Actual driving range will vary based on driving and battery conditions.

\*\*Includes weight of base, curved backrest and batteries.

## Colour Options

