

teal


Roku® Bedside Cabinet

Extremely durable one-piece moulded bedside cabinet to accompany Roku bed. Ideal for diverse environments including high-security and mental health care.



**Shown in Iron Gate with accompanying Roku Bed.*

- Smooth, high impact polyethylene, continuously moulded for maximum durability
- Stain-resistant, seamless, and easy to clean for infection management
- Bedside cabinet sides fit flush to walls
- Smooth, round profile edges ensure enhanced safety and security
- 100% Recyclable material



Discover a spectrum of colours | [Contact us for more information](#)

BEA828500

Roku® Bedside Cabinet

Overall
Height

595mm

Overall
Width

510mm

Overall
Length

460mm

Product
Weight

12.5kg