

The Perfect Powered Wheelchair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.*



When not in use, the Model F can be conveniently stored in a closet or out-of-the-way corner of the room.





800mm

465mm

555mm

855mm

*Model F is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

Lightweight

Weighing only 24kg, the Model F is one of the lightest power chairs on the market today.**







Top Speed

Range

Obstacle Clearance



780mm turning radius for tight turns

Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.







Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



^{**}Weight is 26.7kg, if battery is not removed. Battery weighs 2.7kg.

Note: Driving range of 20kms will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

Customisation

Choose from five colours and a great selection of accessories to match your personal style and enhance your driving experience.

Colour Options ○ White Black WCA870720 (Standard)

Red WCA870730

Controller can be positioned

on either the left or right side.

Liaht Blue

WCA870740

Light Green WCA870750

Accessories

Smart Key WCA870700



Cane Holder WCA870710



Lap Belt WCA870420



Battery* WCA870400



Battery Charger* WCA870410

Specifications

Product Code PWS655500 **Driving Range** Up to 20km Foldable Yes Turning Radius 780mm Seat Sizes Available 450mm / 18in Device Width 555-605mm Device Weight 26.7kg App Enabled Yes, with remote control Max Weight Capacity Up to 113kg Battery Lithium-Ion

Maximum Speed" 8km/h 10° Max Incline Obstacle Clearance 35mm

Width Between Adjustable: 2 positions (400-450mm) Armrests

Device Length 935mm (465 when folded)

~ 5 hours Charging Time

Air Travel Approved Yes Swappable Arm

Covers 5 Colours

Water Resistance IPX4

^{1.} This device is intended to provide mobility to individuals who are comfortable with operating a powered wheelchair. Persons who find it difficult to use a powered wheelchair should refrain from using the device. 2. The device is designed for use on ordinary paved roads or indoor use.



1300 133 120 | aidacare.com.au





(25.3V 10Ah)

^{*}Come as standard

^{***}Acceleration and/or deceleration may change when driving on steep slopes.