

## K Care Shower Commode Weight Bearing Platform



- Assists transportation of the user for showering and can be positioned over the toilet for toileting
- Weight bearing platform footrest and swing out arms assists safe patient transfers in and out of commode
- Ideal for multi user situations such as nursing facilities and hospitals
- Padded or plastic flip up toilet seat available
- Mounted on 4 castors (2 are locking for stability)
- Stainless steel frame for corrosion resistance
- \*Self Propelled model also available
- Note: Pictured with an optional closed padded commode seat (BTC075500)
- **SEATS, PANS & PAN-CARRIERS ARE SOLD AS SEPARATE ITEMS TO THE COMMODE**



### Compliance

AS/NZS 3973:2009

### Warranty Information

- 1 Year

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTC070650	Attendant Propelled Commode - Weight Bearing Platform - 440 mm Wide	430 mm	440 mm	590 mm	150 kg	12 kg
BTC077950	Self Propelled Commode - Weight Bearing Platform - 440 mm Wide	430 mm	440 mm	590 mm	150 kg	17 kg